# TWO CHANNEL CLINICAL AUDIOMETER

## AUDIOSTAR PRO **INSTRUCTIONS FOR USE**

# Table of Contents

EN - English

CS - Czech

FR - French

HU - Hungarian

IT - Italian

TR - Turkish

# TWO CHANNEL CLINICAL AUDIOMETER

## AUDIOSTAR PRO **INSTRUCTIONS FOR USE**

**Title:** GSI AudioStar Pro Instructions for Use

**Manufacturer**


Grason-Stadler  
10395 West 70th Street  
Eden Prairie, MN 55344  
USA

**Copyright © 2021 Grason-Stadler.** All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any means without the prior written permission of Grason-Stadler. The information in this publication is proprietary to Grason-Stadler.

This product may incorporate software developed by others and redistributed according to license agreements. Copyright notices and licenses for these external libraries are provided in a supplementary document that is included with the other software applications distributed with this product.

**Compliance**

The CE 0123 mark identifies compliance with the Medical Device Regulation (EU) 2017/745. Grason-Stadler is an ISO 13485 certified corporation.

European Authorized Representative

Grason-Stadler  
c/o DGS Diagnostics A/S  
Audiometer Alle 1  
5500 Middelfart  
Denmark

**0123**

**Caution:** US Federal law restricts this device to sale by or on the order of a physician or licensed hearing care professional.

---

## TABLE OF CONTENTS

Table of Contents .....	1
Preface.....	7
Manual Conventions.....	7
Regulatory Symbols.....	8
Device Symbols .....	9
Warnings and Cautions.....	11
Important Safety Instructions.....	11
Safety Precautions .....	11
Fluids .....	12
Damage.....	12
Connecting Additional Equipment.....	12
EXPLOSION .....	13
ENVIRONMENTAL FACTORS.....	13
ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS.....	13
Introduction .....	14
Indications for Use.....	14
Intended Use .....	14
Contraindications.....	14
Essential Performance.....	14
Unpacking the system.....	15
Standard Components .....	16
Getting Started.....	17
System Assembly .....	17
Base Unit Rear Panel .....	17
Connections.....	17
Right Side Panel .....	20
Power Switch .....	20
Monitor Speaker.....	20

Left Side Panel.....	20
USB Port.....	21
A/B Cable.....	21
Front Panel .....	22
Monitor .....	22
Elimination of Ambient Noise .....	23
Maximum Ambient Noise (dB re: 20 $\mu$ Pa to the nearest 0.5 dB) .....	23
Operating the AudioStar Pro .....	24
Preliminary Checks .....	24
Prepare to Test.....	24
Placement of the Earphones.....	25
Placement of the Insert Phone .....	25
Placement of the High Frequency Transducer .....	25
Placement of the Bone Vibrator.....	25
Loud Auditory Output.....	25
Test Type Screens .....	26
Title Bar .....	26
Test Type Information.....	26
Navigation Menu .....	27
Time and Date .....	27
Common Icons .....	27
Tone Test Type .....	28
High Frequency Testing.....	29
TEN Test Type.....	29
Presentation Level.....	29
Test Instructions.....	30
Test Procedure .....	30
Interpretation .....	30
Speech Test Type Button.....	30
Integrated Word Files.....	31

More Test Type button .....	32
BKB-SIN .....	33
Presentation Level.....	33
Test Instructions.....	33
Test Procedure .....	34
QuickSIN .....	35
Presentation Level.....	35
Test Instructions.....	35
Test Procedure .....	36
Configuration and Update .....	37
Update.....	37
Bone .....	37
Assistant Headset.....	38
Export Logs .....	38
License.....	38
Settings.....	38
Cleaning and Maintenance .....	39
Routine Maintenance .....	39
Biological Calibration Check .....	39
Periodic Checks .....	39
Earphone and Bone Vibrator Cords .....	39
Hum and Noise.....	39
Distortion and Frequency Shift .....	39
Speech Level Check.....	40
Internal Controls Check .....	40
Bone Vibrator Check .....	40
Masking Level Check .....	40
Talk Forward Check.....	40
Calibration .....	40
Cleaning the System .....	41

Cleaning and Disinfecting Agents .....	41
Single-use Accessories .....	41
Status and Error Messages .....	42
Error Messages .....	42
Please try another selection .....	42
Error .....	42
Error .....	42
Supplies and Accessories .....	43
Appendix A – Technical Data .....	44
Appendix B – Specification of Input/output Connections .....	50
Mains Cable .....	50
Headphones/Inserts (2 connectors, R, L) .....	50
Bone Oscillator .....	50
Gooseneck Mic .....	50
Patient Response x2 .....	51
Talk Forward, Talk Back Mic .....	51
Monitor Headset .....	51
Aux Headset .....	51
External A+B .....	52
VRA .....	52
Monitor .....	52
Speakers (FF1, FF2, FF3, FF4) .....	53
USB A (Host) x4 .....	53
USB B (device) .....	54
Speaker 1 &2 (5 pole DIN) .....	54
LAN Ethernet .....	54
Appendix C - Calibration Reference & Maximum Levels .....	55
Earphones - Pure Tone RETSPL .....	56
Earphones - ANSI Speech RETSPL .....	57
Earphones - IEC Speech RETSPL .....	58

---

Earphones - Pure Tone max HL .....	59
Earphones - NB noise effective masking level.....	60
Earphones - NB noise max HL .....	61
Earphones - ANSI Speech max HL .....	62
Earphones - IEC Speech max HL.....	62
Insert Earphones - Pure Tone RETSPL.....	63
Insert Earphones - ANSI Speech RETSPL.....	64
Insert Earphones - IEC Speech RETSPL .....	64
Insert Earphones - Pure Tone max HL .....	65
Insert Earphones - NB noise effective masking level.....	66
Insert Earphones - NB noise max HL .....	67
Insert Earphones - ANSI Speech max HL .....	68
Insert Earphones - IEC Speech max HL.....	68
Bone Vibrator - Pure Tone RETFL .....	69
Bone Vibrator - ANSI Speech RETSPL .....	70
Bone Vibrator - IEC Speech RETSPL.....	70
Bone Vibrator - Pure Tone max HL.....	71
Bone Vibrator - NB noise effective masking level .....	72
Bone Vibrator - NB noise max HL.....	73
Bone Vibrator - ANSI Speech max HL.....	74
Bone Vibrator - IEC Speech max HL .....	74
Free Field Speakers – ANSI RETSPL and Max HL .....	75
Appendix D - Warranty and Repair .....	76
Appendix E – Recycling and Disposal.....	77
Appendix F – Electromagnetic Compatibility (EMC).....	78
Cautions regarding EMC.....	78
Guidance and manufacturer’s declaration Electromagnetic emissions.....	79
Recommended separation distances between portable and mobile RF communications equipment .....	80
Guidance and Manufacturer’s Declaration Electromagnetic Immunity .....	81

## PREFACE

This user manual provides information about the GSI AudioStar Pro. This manual is intended for technically qualified personnel. This manual contains information pertinent to the use of the AudioStar Pro system including safety information as well as maintenance and cleaning recommendations. Proper use of this device depends on careful reading of all instructions and labeling. Also be sure to follow all safety standards set by your facility.

**Please note:** This User Manual is not intended as a training manual. The reader should consult standard audiology texts for the theory and application of the tests provided by this instrument.

READ THIS ENTIRE MANUAL BEFORE ATTEMPTING TO USE THIS SYSTEM!

## MANUAL CONVENTIONS

Throughout this manual, the following meaning of warnings, cautions and notices are used.

### WARNING

The WARNING symbol identifies conditions or practices that may present danger to the patient and/or user.

### CAUTION

The CAUTION Symbol identifies conditions or practices that could result in damage to the equipment

---

**NOTE:** Notes help you identify areas of possible confusion and avoid potential problems during system operation.

---

REGULATORY SYMBOLS

Symbol	Description
	Conforms to European Medical Device Regulation (EU) 2017/745.
	Symbol for "SERIAL NUMBER."
	GSI Part Number.
	Return to Authorized Representative, Special disposal required.
	Medical Equipment Classified by Intertek Testing Services NA Inc. with respect to electric shock, fire, and mechanical hazards only, in accordance with UL 60601-1. Classified under the Medical Device Regulation (EU) 2017/745 as a Class IIa device.
	China RoHS symbol for product compliance
	Symbol for "European Authorized Representative."
	Symbol for "Manufacturer."
	Symbol for "Date of Manufacture."
	On/Off
	Keep Dry

Symbol	Description
	<p>This side up.</p>
	<p>Consult the operating instructions/directions for use. A copy of the operating manual is available on this website: <a href="http://www.grason-stadler.com">www.grason-stadler.com</a> A printed copy of the operating instructions can be ordered from Grason-Stadler for shipment within 7 days; or you can contact your local representative.</p>
	<p>Consult the operating instructions/directions for use. A copy of the operating manual is available on this website: <a href="http://www.grason-stadler.com">www.grason-stadler.com</a> A printed copy of the operating instructions can be ordered from Grason-Stadler for shipment within 7 days; or you can contact your local representative.</p>

## DEVICE SYMBOLS

The following symbols appear on the instrument.

**Definition:** Type B applied part – IEC 60601-1 uses the term applied part to refer to the part of the medical device which comes into physical contact with the patient for the device to carry out its intended function. This icon indicates that the GSI AudioStar Pro is Type B according to the international standard IEC 60601-1. The applied parts are the headphones, insert earphones, and the bone oscillator.

Symbol	Description
	Consult Operating Instructions.
	Follow Instructions for Use.
	Connection for headphones. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Connection for the insert earphones. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Connection for the bone conductor transducer.
	Connection for the sound field speakers. The red markings around the connector indicate the right ear and the blue marking indicate the left connector.
	Monitor
	Patient response switch
	On/Off - Next to power mains.

## WARNINGS AND CAUTIONS

### IMPORTANT SAFETY INSTRUCTIONS

The following safety precautions must always be observed. General safety precautions must be followed when operating electrical equipment. Failure to observe these precautions could result in damage to the equipment and injury to the operator or patient.

The instrument must only be used by hearing health care professionals qualified to perform hearing threshold and auditory tests such as an audiologist, otolaryngologist, researcher, or a technician under direct supervision by the specialist. Users should use their professional skills when interpreting the results and this should be done in conjunction with other testing as deemed appropriate given their professional skills. Incorrect use could lead to wrong results.

The employer should instruct each employee in the recognition and avoidance of unsafe conditions and the regulations applicable to his or her work environment to control or eliminate any hazards or other exposure to illness or injury. It is understood that safety rules within individual organizations vary. If a conflict exists between the material contained in this manual and the rules of the organization using this instrument, the more stringent rules should take precedence.

***In case of death or serious incident in relation to the use of the device, the incident must immediately be reported to Grason-Stadler and the local national competent authority.***

### SAFETY PRECAUTIONS

This product and its components will perform reliably only when operated and maintained in accordance with the instructions contained in this manual, accompanying labels, and/or inserts. A defective product should not be used.

Make sure all connections to external accessories are snug and secured properly. Parts which may be broken or missing or are visibly worn, distorted, or contaminated should be replaced immediately with clean, genuine replacement parts manufactured by or available from GSI.

This instrument is suitable in hospital environments except for near active high frequency (HF) surgical equipment and radio frequency (RF) shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

In case of emergency, disconnect the instrument from the supply mains by pulling the plug out of the power socket.

Before the first use of the instrument each day, or if suspect or inconsistent results are apparent, the checks specified in the Routine Maintenance section should be carried out. If the system is not functioning properly, do not operate it until all necessary repairs are made and the unit is tested and calibrated for proper functioning in accordance with Grason-Stadler published specifications.

No parts of the equipment can be serviced or maintained while in use with the patient.

Latex is not used anywhere in the manufacturing process. The base material for the earphone cushions is made from natural and synthetic rubber.

### Fluids

The device is not intended to be used in environments exposed to fluid spills. No means is specified for fluid protection (not IP classed). Do not use the device in the presence of fluid that can contact any of the electronic components or wiring. Should the user suspect fluids have contacted the system components or accessories, the unit should not be used until deemed safe by an authorized service technician. Do not immerse the unit in any fluids. See the Routine Maintenance Section of this manual for the proper cleaning procedure for the instrument and its accessories and the function of single-use parts.

### Damage

Do not drop or otherwise impact this instrument. If the instrument is dropped or damaged, return it to the manufacturer for repair and/or calibration. Do not use the instrument if any damage is suspected.

Do not attempt to open, modify or service the instrument. Return the instrument to the manufacturer or distributor for all repair and servicing requirements. Opening the instrument will void the warranty. See Warranty and Repair in the Appendix

### Connecting Additional Equipment

This equipment is intended to be connected to other equipment thus forming a Medical Electrical System. External equipment intended for connection to signal input, signal output or other connectors shall comply with the relevant product standard e.g., IEC 60950-1 for IT equipment and the IEC 60601-series for medical electrical equipment. In addition, all such combinations – Medical Electrical Systems – shall comply with the safety requirements stated in the general standard IEC 60601-1, edition 3, clause 16. Any equipment not complying with the leakage current requirements in IEC 60601-1 shall be kept outside the patient environment i.e., at least 1.5 m from the patient support or shall be supplied via a separation transformer to reduce the leakage currents.

Any person who connects external equipment to signal input, signal output or other connectors has formed a Medical Electrical System and is therefore responsible for the system to comply with the requirements. If in doubt, contact a qualified medical technician or your local representative. If the instrument is connected to a PC (IT equipment forming a system) ensure not to touch the patient while operating the PC. If the instrument is connected to a PC (IT equipment forming a system) assembly and modifications shall be evaluated by qualified medical technician according to safety regulations in IEC 60601.

## EXPLOSION

This system is not explosion proof.

Do NOT use in the presence of flammable gaseous mixtures. Users should consider the possibility of explosions or fire when using this device near flammable anesthetic gases.

Do NOT use the system in a highly oxygen-enriched environment, such as a hyperbaric chamber, oxygen tent, etc.

## ENVIRONMENTAL FACTORS

Use and store the instrument indoors only. It is recommended that the instrument be operated within an ambient temperature range of 15 °C / 59 °F to 35 °C / 95 °F and in relative humidity between 30 % and 90 % (non-condensing).

Transport the instrument in temperatures between 20 °C / -4 °F to +50 °C / +122 °F and store the instrument in temperatures between 0 °C / 32 °F to 50 °C / 122 °F.

## ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information in the Appendix. Although the instrument fulfills the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g., from mobile phones etc. This appendix provides guidance on the electromagnetic environment in which to operate the instrument.

Portable and mobile radio frequency (RF) communications equipment can affect medical electrical equipment. Use of this instrument adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

## INTRODUCTION

Thank you for purchasing the GSI AudioStar Pro™ system. The GSI AudioStar Pro is a stand-alone, two channel, advanced clinical/research audiometer. It supports data transfer to a computer but does not rely on the computer to perform tests and will not lose data if the computer connection fails. User login and password controls provide security for patient data. Test data can be transferred to software such as GSI Suite and Noah and can be integrated into your EMR/EHR system.

## INDICATIONS FOR USE

The AudioStar Pro is intended to be used for the identification of hearing loss in patients of any age. It is intended to be used by an audiologist, ENT, hearing healthcare professional, or trained technician in a hospital, clinic, healthcare facility or other suitable quiet environment as defined in ANSI S3.1 or equivalent.

## INTENDED USE

The GSI AudioStar Pro is intended to be used for the purposes of determining patient hearing sensitivity. It is intended to quantify the patient's level of hearing by presenting pure tones stimuli through specific transducers at different frequencies and at different sound pressure levels.

## CONTRAINDICATIONS

There are no known contraindications for use.

## ESSENTIAL PERFORMANCE

Essential Performance for this instrument as per IEC 60601-1 is defined by the manufacturer as:

This instrument does not have an Essential Performance.

Absence or loss of Essential Performance cannot lead to any unacceptable immediate risk. Final diagnosis shall always be based on clinical knowledge.

## UNPACKING THE SYSTEM

Although this GSI AudioStar Pro Clinical Audiometer was carefully tested, inspected, and packed for shipping, it is good practice after receiving the instrument to immediately examine the outside of the container for any signs of damage. Notify the carrier if any damage is observed.

- ✓ It is recommended that you unpack your AudioStar Pro carefully to make sure that all components are removed from the packing materials.
- ✓ Verify that all components are included as shown on the packing slip included with your shipment.
- ✓ If any component is missing, contact your distributor immediately to report the shortage.
- ✓ If any component appears to be damaged in shipment, contact your distributor immediately to report it. Do not attempt to use any component or device that appears to be damaged.
- ✓ Please check that all accessories listed below have been received in good condition.
- ✓ All standard components are part of the medical device and suitable for operation in a patient environment.

Notify the carrier immediately if any mechanical damage is noted. This will ensure that a proper claim is made. Save all packing material so the claim adjuster can inspect it as well. Notify your dealer or GSI when the adjuster has completed the inspection.

***Save all the original packing material and the shipping container so the instrument can be properly packed if it needs to be returned for service or calibration.***

## STANDARD COMPONENTS

- AudioStar Pro Base Unit and power cord
- Subject Response Hand Switch
- DD45 Audiometric Headset P3045 (Adult) for DD45 systems
- DD65v2 Audiometric Headset for DD65v2 systems
- DD450 Audiometric HF Headset for DD450 and AMTAS systems
- IP30 Insert earphones, 10 ohm
- Bone Vibrator Assembly, B81, 10 ohm
- Monitor Headset with boom mic
- Headphones, Assistant (Aux Intercom) with inline control
- Extension cable for Assistant headset (stereo, m-f, 3.5mm)
- Gooseneck Microphone
- Talk Back Microphone (TB-2) with mounting bracket
- Wireless keyboard and mouse
- Earpads for monitor headset with boom mic
- Microphone windscreen (Sennheiser PS 01 for PC 131 headset)
- Cable, USB A/B, 2 meters
- Patch Cord, 6 feet, red, 1 each
- Patch Cord, 6 feet, blue, 1 each
- Patch Cord, 6 feet, black, 1 each
- Dust Cover
- Remote Keyboard Shortcut Guide
- AudioStar Pro Software and Manuals, USB Bundle
- GSI Suite Data Management Software, USB Bundle

## GETTING STARTED

### SYSTEM ASSEMBLY


1. Connect transducers to audiometer base unit
2. Connect USB cable from base unit to PC if using GSI Suite
3. Connect power cable from base unit to outlet or isolation transformer

### BASE UNIT REAR PANEL


#### Connections

The connectors on the rear panel of the GSI AudioStar Pro are shown below. The label and jacks are visible by turning the instrument around on a flat, stable surface.

	Connection	Description	Graphic
1	<b>Mains Power Input</b>	IEC 14	
2	<b>USB Computer Connection</b>	USB B style connector	
3	<b>LAN Connections</b>	Ethernet Connection RJ45 Currently not supported	
4	<b>USB Connectors</b>	USB A style plug	

	Connection	Description	Graphic
5	<b>External Monitor Output</b>	HDMI Video only signals, no audio, 1024 x 768 resolution	
6	<b>FF Speaker DIN Connection Output</b>  SFS - Sound Field Speaker	5 pin DIN connector Provides connection between the internal amplifier to left and right loudspeakers within a sound room  <b>NOTE:</b> Free Field Speaker Outputs 1 and 2 are 25 Watts per channel into 8 ohm.	
7	<b>FF Speaker RCA Connections Output</b>	4 RCA jacks Optionally connect to 4 speakers through an external amplifier using jacks 1 -4 (contact a GSI Representative for more information)  <b>NOTE:</b> Free Field Line Outputs 1 and 2 are 5 VRMS into a 2000 ohm load.  <b>NOTE:</b> Cannot use internally amplified speaker connection and externally amplified speaker connections at the same time.  <b>NOTE:</b> The CTL connection is for future use - not currently supported.	
8	<b>Left and Right Headphone Outputs</b>	6.35 mm stereo jack Left (blue) and Right (red)	
9	<b>Patient Response Inputs</b>	6.35 mm mono jack 1 or 2 hand switches can be used	
10	<b>AUX Intercom Output</b>	3.5mm stereo jack Assistant monitor headset connector	

	Connection	Description	Graphic
11	<b>Left and Right Insert phone Outputs</b>	6.35 mm stereo jack Left (blue) and Right (red)	
12	<b>VRA Connection Output</b>	3.5 mm TRRS jack to activate a left, center, or right VRA system. (Ref. Appendix for pinouts) (contact a GSI service representative for details)	
13	<b>Left and Right High Frequency Headset Output</b>	6.35 mm stereo jack Left (blue) and Right (red)	
14	<b>Talkback Microphone Input</b>	3.5 mm stereo jack <b>NOTE:</b> Microphone inputs are between .25 mV and 5 mV for a 0 dB reading on a VU indicator; the input impedance is 3,200 ohm.	
15	<b>Ext. A and B</b>	3.5 mm stereo jack Input jacks for optional digital music player or CD player input <b>NOTE:</b> External A and B inputs are between 15 mV and 500 mV for a 0 dB reading on a VU indicator; the input impedance is 50,000 ohm.	
16	<b>Bone Vibrator</b>	6.35 mm phone stereo jack	

## RIGHT SIDE PANEL

### Power Switch

The power switch is located on the right-side panel. To turn on/off the GSI AudioStar Pro, use the power switch.

**NOTE:** Do not block access to the power switch.

### Monitor Speaker

The monitor speaker is located on the right-side panel. The monitor speaker can be activated by pressing the monitor speaker button on the front panel. The level of the Channel 1 and Channel 2 stimuli can be adjusted using the monitor knob on the front panel of the instrument.



## LEFT SIDE PANEL

The following connectors are found on the left side panel of the GSI AudioStar Pro.

Connection	Description	Graphic
USB Ports	2 USB ports (A style)	

Connection	Description	Graphic
<b>Monitor Microphone</b>	3.5 mm stereo jack Monitor microphone	
<b>Headphones</b>	3.5 mm stereo jack Monitor headphones	
<b>Gooseneck Microphone</b>	6.35 mm stereo jack (optional)	

### USB Port

The AudioStar Pro is equipped with four (4) USB ports. It is possible to connect external devices such as mouse, keyboard, or external printer to be used with the audiometer. Additionally, a memory stick can be inserted into a USB port for updating software, adding additional sound files, license updates, printing to PDF or exporting diagnostic log files.

---

**NOTE:** Scan files on a USB drive for viruses prior to installing the drive into the instrument.

---

### A/B Cable

Remote connection to an external computer is achieved using a standard A/B USB cable.


---

**NOTE:** It is recommended to always have the USB ports enabled on the PC. Disable the "suspend USB" option on the PC.

---

## FRONT PANEL

The front panel controls are used to operate the device. The blue shaded areas have duplicated controls that operate the 2 channels. The areas outside the blue shaded area control the functions of the device that are channel independent. The controls are labeled to indicate their function. The four navigation buttons, and the middle select button, can be used to make selections from the on-screen menus as well as navigate through the internal files for speech testing.

## MONITOR

The AudioStar Pro comes standard with an LCD display. The LCD is hinged to the GSI AudioStar Pro and is used to display all the testing information from the instrument. The information displayed on the AudioStar Pro LCD varies depending on the Test Type.

## ELIMINATION OF AMBIENT NOISE

Excessive noise in the test environment reduces test validity because it masks the test signals. Maximum permissible ambient noise levels (MPANLs) are defined in the standard ANSI S3.1-1991 (R2008) American National Standard Maximum Permissible Ambient Noise Levels for Audiometric Test Rooms.

The following table shows the maximum background levels that can be present inside the room while a valid hearing test from 125 Hz to 8000 Hz is being conducted. These values apply for hearing threshold measurements to 0 dB HL. The standard specifies that MPANLs should be adjusted appropriately when hearing thresholds for pure tones are measured above and below 0 dB HL.

Maximum Ambient Noise (dB re: 20  $\mu$ Pa to the nearest 0.5 dB)

	Supra-aural Earphone	Insert Earphone	Ears not covered i.e., Sound Field
<b>125</b>	35	59	29
<b>250</b>	25	53	21
<b>500</b>	21	50	16
<b>1000</b>	26	47	13
<b>2000</b>	34	49	14
<b>4000</b>	37	50	11
<b>8000</b>	37	50	14

**NOTE:** Live voice testing requires a separate sound attenuated room for the patient to avoid feedback and direct transmission of the test stimuli.

Sound Attenuation (dB) for Earphones per ISO 4869-1

	DD45 Earphone	Insert Earphone	DD450
<b>125</b>	3	33.5	14.5
<b>250</b>	5	34.5	16.0
<b>500</b>	7	34.5	22.5
<b>1000</b>	15	35.0	28.5
<b>2000</b>	26	33.0	32.0
<b>4000</b>	32	39.5	45.5
<b>8000</b>	24	43.5	44.0

## OPERATING THE AUDIOSTAR PRO

### PRELIMINARY CHECKS

Before starting any procedures using the GSI AudioStar Pro Clinical Audiometer, ensure that the power cord is plugged into a properly grounded receptacle.

Check also that all cords from the transducers, patient response hand switch (if used), and printer fit securely in their connectors on the rear and side panels. Inspect all cords for fraying and damage. If there is any damage to any cord, do not use the AudioStar Pro. If speech testing with recorded voice from an external source is to be performed, check that the CD or digital music player device is connected and operating properly.

It is recommended that all parts that come into direct contact with the patient (e.g., earphone cushions) are subjected to standard disinfecting procedures between patients. This includes physically cleaning and using a recognized disinfectant. Individual manufacturer's instructions should be followed for use of any disinfecting agent to provide an appropriate level of sterilization.

Handle earphones, bone vibrator, and insert earphones with care. Do not drop them nor allow them to be banged together. Severe mechanical shock can alter their operating characteristics or change the output levels, which may require that the transducers be replaced.

### PREPARE TO TEST

Once you have performed the preliminary checks and any routine maintenance checks prepare to test patients.

- √ Turn on the instrument and allow it to come to operating temperature (approximately 10 minutes).
- √ Check that the transducers and other system components are operating properly.
- √ Seat the patient comfortably in the test area.
- √ Place the selected transducers on the patient.

Do not connect or disconnect Earphones, Insert Phones, Bone Vibrator, High Frequency transducers, or any other accessories while in contact with the patient.

### Placement of the Earphones

Prior to positioning the earphones on the patient's head, inspect the ear canals for any blockage due to cerumen or foreign objects. Recognize that soft-walled ear canals may collapse under the earphones, and this can lead to incorrect threshold levels. Insert phones might be used in these cases. Eliminate all obstructions, such as glasses, hair, or hearing aid, between the earphone and the patient.

Center the earphones over both ears and adjust the headband so that it rests solidly on the crown of the head and exerts pressure on both ears. Place the earphone with the red connector over the patient's right ear and the earphone with the blue connector over the left ear.

### Placement of the Insert Phone

Be sure there is an eartip attached to the insert phone before inserting into the patient's ear. Inserting the insert phone without an eartip could cause harm to the patient. Push the correctly sized eartip onto the earphone and then place the insert phone securely into the patient's ear. When using the insert phones, follow the manufacturer's recommended procedure for eartip preparation, placement, and insertion.

Insert eartips are single use only. Using disposable eartips ensures sanitary conditions for each patient. Failure to comply with the single-use requirement can lead to cross-contamination and infection.

### Placement of the High Frequency Transducer

Remove eyeglasses and earrings if possible and position the transducer directly on the head of the patient. Place the rubber cushions so that the earphone diaphragm is aimed directly at the opening into the ear canal. Adjust the headband for a tight fit. If the cushions are not tight to the ears, the test result will be false, especially at lower frequencies.

### Placement of the Bone Vibrator

The bone vibrator can be placed on the promontory of the mastoid process or on the forehead, whichever has been selected in the configuration application or modified in the Configure screen.


### Loud Auditory Output

The AudioStar Pro audiometer is capable of outputting loud sounds and is designed to provide a visual indicator to the examiner when presenting sounds at or above 100 dB. The level display is shown with a yellow background when loud output is selected for presentation.

## TEST TYPE SCREENS

The information displayed on the AudioStar Pro LCD varies depending on the Test Type. There are common elements found on all screens such as the Channel 1 and 2 level settings, the Navigation menu, and the Title Bar.

### Title Bar

The title bar is located at the top of the display. The title bar displays the test type in the middle. The patient's name will appear on the left side of the title bar if a patient name has been entered (or selected from an imported patient list). The right side of the title bar displays the examiner name if examiners have been entered. The examiners can be entered from the Configuration application.

### Test Type Information

Under the title bar test specific information will be displayed. On the left and right side, the current output in dB HL for Channel 1 and Channel 2 will be displayed. The other information displayed will depend on the test type and is described as part of the individual test type displays.

## Navigation Menu

This menu is located at the bottom of the display. It utilizes the on-board navigation buttons or an external mouse to access the menu options. The menu is specific to the test type selected.

## Time and Date

The date and time are displayed in the bottom right corner of the screen. Using the Configuration Application, the Time can be configured in a 12- or 24-hour format and the Date can be configured in any order (dd/mm/yyyy, etc.). It is also possible to update the date and time in the configuration screen of the instrument. Use the navigation buttons and the attenuator dials to set the appropriate time and date. The date and time will be updated when a new session is started.

---

**NOTE:** The time does not change automatically for daylight savings time. The operator must manually change the time using the configure button on the front panel of the instrument or the configuration application.

---

## Common Icons

These icons are found in the test information area and common to the different test types.

**Talk Forward** – When pressed, a head with a headset icon will appear. This icon will remain active if the talk forward button is depressed.

**Store** – When either of the store buttons is pressed, a floppy disc icon flashes and the result is then displayed.

**Interlock** – When interlock is active, a padlock icon will appear.


**Tracking** – When tracking is selected, a railroad track icon will appear.

**Aux Intercom** – When pressed, the Aux intercom icon indicates direct communication between the operator and the Aux headset.

**Data Transfer** – When there is an active connection between the AudioStar Pro and an external computer, communication will be indicated by the blue arrows.

**Left, Center and Right VRA** – A VRA system with a right, left and centering toy can be connected to the AudioStar Pro. The VRA-L, VRA-R and VRA-C icons will appear on the display when the front panel buttons have been pressed to activate the VRA system.


**Pencil Icon** - This icon opens a comments window (must use external keyboard to utilize comment section). Comments can be entered from any test screen, and it is possible to review and edit comments from any test screen.

## TONE TEST TYPE

Pressing the Tone Test Type button prepares the GSI AudioStar Pro™ for pure tone air and bone conduction testing from 125 – 20,000 Hz.

Use the Range menu to select the desired frequency range to be tested: Standard (125-8,000 Hz), High (8,000- 20,000 Hz), or Full (125 – 20,000 Hz). Each selection on the blue navigation menu is specific to Pure Tone Testing. Pressing this button will set the defaults from the configuration application to start the test type. Auto Hz can be defined in the Config App and automatically moves the next designated test frequency when “store” is pressed.

---

**NOTE:** If you select a different transducer or stimulus or ear and leave the tone test type, when you return to tone, the AudioStar Pro will default to the last settings that were selected.

---

- Press the Tone Test Type Button.
- Verify that the transducers and signals are correct.
- Verify the frequency range is correct.
- Perform air conduction threshold testing.
- Press “store” after each threshold is obtained.

When the pure tone evaluation is complete, move to the next test type in your typical testing sequence.

## HIGH FREQUENCY TESTING

Use the navigation buttons or the mouse to access the Range menu. Available frequency ranges are Standard (125-8,000 Hz), High (8,000-20,000 Hz) or Full (125 – 20,000 Hz). Auto Hz can be defined in the Config App and automatically moves the next designated test frequency when “store” is pressed.

- Press the Tone Test Type Button.
- Select the desired frequency range (High or Full).
- Verify that the transducers and signals are correct.
- Perform High Frequency Testing.
- Press “store” after each threshold is obtained.

When the high frequency evaluation is complete, move to the next test type in your typical testing sequence

## TEN TEST TYPE

### Presentation Level

- If the audiometric threshold is 60 dB or less, set the TEN level to 70 dB.
- If the audiometric threshold is 65 dB or greater, set the TEN level 10 dB higher than the threshold.
- If the TEN noise is reported to be too loud, set the TEN level at 90 dB (maximum) or the same level as audiometric threshold.

### Test Instructions

When the starting level has been determined, instruct the patient in the same manner as when measuring pure tone thresholds with masking.

### Test Procedure

Using the blue navigation buttons, access the TEN test type to prepare the AudioStar Pro for TEN testing. The procedure for determining thresholds in the TEN is identical to the manual pure tone audiometry except that a 2 dB final step size should be used for maximum accuracy. The TEN will take approximately 4 minutes per ear (to complete all test frequencies).

### Interpretation

The accepted rule is that a dead region is present when the TEN-masked threshold is at least 10 dB above the audiometric threshold and the TEN Threshold is 10 dB above the TEN noise.

## SPEECH TEST TYPE BUTTON

Pressing the Speech Test type button prepares the AudioStar Pro for Speech testing. The internal .Wav files can be presented by either using the present button or by a single click of a wireless mouse. The correct/incorrect/clear buttons can be used to score. It is also possible to automatically play wordlists by pressing the interrupt button. These preferences are determined in the Config App. It is critical that the test type be carefully selected as the reporting/storing is dependent upon test type. If you wish to do a PIPB rollover, you can select the speech audiogram page.

## Integrated Word Files

When Speech Test Type is selected, the AudioStar Pro defaults to internal speech files. These can be presented for consistent recorded speech testing.

- Define word list favorites in the Config App for fast selection of common lists.
- Utilize the navigation menu or external mouse to select the test type and the wordlist.
- Select Word Nav and use the navigation buttons to highlight word stimulus. Press the present bar to present the word.


-OR-

- Utilize the navigation buttons external mouse to present the words (single click to present).
- When the speech stimulus is being presented, the word will be highlighted yellow.

- When the patient responds (and the yellow highlight disappears), the stimulus word/sentence can be scored correct or incorrect.
- The stimulus word/sentence will turn green for correct or orange/red for incorrect. The center area of the display will indicate the % correct/#words presented.
- After the completion of each speech test type, press store to save the results in the speech results table.

When your speech evaluation is complete, move to the next test type in your typical test sequence.

### MORE TEST TYPE BUTTON

Pressing the "More" test type button calls up a menu of the following special tests: ABLB, BKB-SIN, QuickSIN, SISI, and Tone Decay. Use the on-board navigation buttons or an external mouse to select the special test.

## BKB-SIN

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' indicators. There are also fields for 'INT A', 'INT B', 'PTA AC: BC:', and 'SII:'. The 'Test Results Group 1' section contains three empty tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. Below this is a 'List Pair 1' section with a table of test sentences and their corresponding scores. The bottom of the interface features a control bar with 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp of 10:03 AM on 02-26-2019.

Score	S/N	Word List	Score
-	21	A1 THEY are LOOKING AT the CLOCK	-
-	18	A2 The CAR ENGINE is RUNNING	-
-	15	A3 CHILDREN LIKE STRAWBERRIES	-
-	12	A4 THEY are BUYING some BREAD	-
-	9	A5 The GREEN TOMATOES are SMALL	-
-	6	A6 HE PLAYED with his TRAIN	-
-	3	A7 The BAG FELL to the GROUND	-
-	0	A8 The BOY DID a HANDSTAND	-
-	-3	A9 The WATER BOILED QUICKLY	-
-	-6	A10 The MAN is PAINTING a SIGN	-
-	Sum		-

### Presentation Level

For standard SNR Loss testing the BKB-SIN Test should be presented at a relatively high level (loud, but below discomfort). Normative data on normal-hearing adults and normal-hearing children were collected using binaural presentation via insert earphones, at a presentation level of 70 dB HL (83 dB SPL). Normative data on adult cochlear implant users were collected using a 65 dB SPL presentation level in sound field (equivalent to 50 dB HL at 0 degrees azimuth).

### Test Instructions

“You will hear a man talking to you through the earphones (or loudspeaker). He is going to say “Ready” and then he’ll say a sentence. Repeat the sentence the man says. You will hear other talkers in the background. Don’t pay any attention to them; just repeat what the man says. The background talkers will get louder, and then it will be hard for you to hear the

man's voice. When that happens, it is OK to guess; repeat anything you think you heard the man say."

### Test Procedure

- Select BKB-SIN from the More Tests Menu.
- Select the proper transducer and intensity levels for each channel.
- Select the appropriate age from the Navigation Menu
- Select the appropriate word list pair. Lists 1-8 are appropriate for all patients. List pairs 9-16 are appropriate for CI patients or those with significant hearing loss.
- Using the Word Nav and front panel navigation buttons or an external mouse, select the first sentence.
- Press the present bar or click the first sentence.
- Score the four/three key words highlighted in each sentence by pressing the INCORRECT button for each word repeated incorrectly by the patient.
- Complete the entire list pair (20 sentences) and press "store".

To interpret the SNR loss score for adults, see table below.

Refer to the BKB-SIN user manual for how to interpret test results for children.

<b>SNR Loss</b>	<b>Degree of SNR Loss</b>	<b>Expected Improvement with Directional Mic</b>
<b>0-3 dB</b>	Normal/near normal	May hear better than normal hearing individuals hear in noise
<b>3-7 dB</b>	Mild SNR loss	May hear almost as well as normal hearing individuals hear in noise
<b>7-15 dB</b>	Moderate SNR loss	Directional microphones help; consider array mic
<b>&gt;15 dB</b>	Severe SNR loss	Maximum SNR improvement is needed; consider FM system

## QUICKSIN

The screenshot displays the Quicksin software interface. At the top, it shows 'Patient ...' and 'Quicksin'. Below this are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 is set to 70 dB HL and Channel 2 to 0 dB HL. Each channel has a frequency scale from -20 to 3 dB HL. Between the channels are 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' tables, each with columns for R, B, and L, and rows for Basic, HFE, and HFE-LP. Below these are 'PTA AC: BC:' and 'SII:' fields for both Right and Left ears, along with a 'Reliability' field. The 'Test Results Group 1' section contains three empty tables with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. At the bottom, there is a 'QuickSin : List 1 (Track 3)' section with a list of sentences and a 'Score' table with columns for S/N and a dash. The 'Score' table shows values for S/N 25, 20, 15, 10, 5, 0, and Sum, all with dashes. The bottom status bar includes 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a timestamp of 10:53 AM 03-01-2019.

### Presentation Level

For pure-tone average (PTA) less than or equal to 45 dB HL, set the attenuator (s) in Channel 1 and Channel 2 to 70 dB HL. For PTA of 50 dB HL or greater, set the attenuators to a level that is judged to be “loud, but okay.” The sound should be perceived as loud, but not uncomfortably loud. Tests can be performed in the right, left or binaural conditions.

### Test Instructions

“Imagine that you are at a party. There will be a woman talking and several other talkers in the background. The woman’s voice is easy to hear at first, because her voice is louder than the others. Repeat each sentence the woman says. The background talkers will gradually become louder, making it difficult to understand the woman’s voice, but please guess and repeat as much of each sentence as possible.”

## Test Procedure

- Select QuickSIN from the More Tests Menu.
- Ensure the proper transducer and intensity levels are selected.
- Select the appropriate word list. Lists 1-12 are standard lists used to determine SNR loss.
- Using the Word Nav and front panel navigation buttons or an external mouse, select the first sentence.
- Press the present bar or click the selected sentence.
- Score the five key words highlighted in each sentence by pressing the INCORRECT button for each word repeated incorrectly by the patient.
- Press Store
- The SNR Loss score will appear in the SCORE/WORD window.

---

**NOTE:** It is recommended that at least 2 lists be presented in each condition. The average score will appear in the Group 1 SNR averages window.

---

To interpret the SNR loss score, see the table above in the BKB-SIN section.

Refer to the QuickSIN user manual for advanced testing options.

## CONFIGURATION AND UPDATE

The configure button presents a Configuration screen where it is possible to view the instrument information such as serial number, software version and last calibration date. This screen displays setup options to update the AudioStar Pro software, configure bone conduction symbol settings, enable/disable the Aux Intercom, export log files, set the date and time, and adjust the screen brightness.

### Update

Place a USB drive with the appropriate update loaded into one of the four USB ports. Select Update and then select from device, settings, or sound files to update the instrument. Software and sound files updates must be obtained from GSI or an authorized GSI representative. Selecting Settings allows for installing an update from a config file generated from the Configuration Application.

### Bone

Select the symbol scheme for bone conduction testing. Choose between MAS (Mastoid) and FOR (Forehead.) This selection will be active throughout the current session. When a new session is started, the symbol scheme will revert to the configured preference.

### Assistant Headset

When the box is checked the Aux Intercom is enabled. If the box is not checked the Aux Intercom is disabled.

### Export Logs

Select to obtain and send the instrument logs to an attached USB drive. The instrument logs may be requested to help troubleshoot any problems.

### License

Select the License icon to see the licensed options on the instrument. In the dialog that is presented, the current options are indicated and there are buttons to allow the update or export of the license key. An update to the license key can be made manually by typing the key code into the text box or the key can also be imported from an attached USB drive.

---

**NOTE:** One license key contains all the instrument options. When a successful update occurs, the option to add Another Key refers to licensed word lists, if applicable.

---

### Settings

Select the settings icon to configure the items listed below.

*Graph Orientation* - Select from the dropdown to control how the audiogram is presented on the screen. Choose left-right, right-left or combined.

*Masked and Unmasked* - Select from the drop-down menu to have masked thresholds replace unmasked thresholds at the same frequency OR store both masked/unmasked at each frequency.

*Printer Protocol* - Select from the dropdown the default printer or printer protocol to be used for the printouts from the instrument. Consult the printer documentation to determine the correct protocol.

*Date and Time* - Select to change the date and update the time displayed on the AudioStar Pro. Use the on-board navigation buttons to tab to each field and the attenuator dials to change the date and time from the stand-alone instrument.

*Brightness* - Select to change the brightness of the screen.

*Save* - Select to save all settings changes.

*Cancel* - Select to cancel all changes to settings.

## CLEANING AND MAINTENANCE

### ROUTINE MAINTENANCE

#### Biological Calibration Check

The design of the GSI AudioStar Pro audiometer should provide trouble-free service for a long time. It is recommended to routinely record and file the audiogram of one person for the purpose of biologic calibration. This person (or group of persons) should have a known stable audiometric curve that does not exceed 25 dB HL at any frequency. This procedure should start when the GSI AudioStar Pro is first installed and then be continued.

Remember that individual thresholds can shift by as much as 5 dB from day to day; however, variations that exceed this range can point to difficulties which require attention.

#### Periodic Checks

The routine maintenance checks described below can detect possible instrument problems. If any issues are found, the instrument should receive technical service before further use. The checks should be made at periodic intervals, even if biologic checks reveal no problems.

#### Earphone and Bone Vibrator Cords

With extended use, all transducer cords tend to fray internally at the connectors. To evaluate the cord status, turn on the GSI AudioStar Pro. Set the HL to a comfortably audible level. Place the transducer on your head. Activate both Interrupt buttons. Bend the cord next to the plug at both ends of each earphone. Listen for an intermittent signal, abrupt changes in the signal level, or a scratchy sound that coincides with the flexing of the cord. The presence of any of these conditions signifies that the cord should be replaced. Repeat this check for all transducers.

#### Hum and Noise

Set the GSI AudioStar Pro to Tone test type with the standard earphones selected and the Channel 1 Interrupt button in the ON mode. Turn the Channel 1 Hearing Level control from 0 to 60 dB HL. Listen for low frequency hum (60 or 120 Hz) and any other noise (hiss or low rushing sound) at all attenuator levels through the earphone. Some audible noise at levels above 70 dB is permissible. If these noises are detected below 70 dB, the audiometer should be scheduled for maintenance. Repeat for Channel 2.

#### Distortion and Frequency Shift

Check for distortion and frequency shift by listening to the GSI AudioStar Pro's output through the earphones at each frequency (in the 125 Hz to 12,000 Hz range) at a loud, but not uncomfortable level (70 to 80 dB HL for normal ears). Listen also to ensure that the signal frequencies change appropriately when the Frequency up arrow (>) and down arrow

(<) pushbuttons are operated. If distortion is heard in one earphone but not the other, the chances are high that the earphones are at fault and should be replaced. In any case, the audiometer should be scheduled for immediate maintenance.

### Speech Level Check

To check the speech level with recorded speech, select the Speech test type button. Place the earphones on a person with normal hearing and present a word list at 40 dB. If intelligible speech is not heard, with the Channel 1 Hearing Level control set at 40 dB, the audiometer should be scheduled for technical service.

### Internal Controls Check

Should the front panel controls lock into one state, and it is not possible to change any of the parameters, turn off the power. Wait one minute and then power on.

### Bone Vibrator Check

This check must be performed in a quiet environment or in a sound room. With the frequency set to 2,000 Hz, the Channel 1 level set at 40 dB HL and the bone vibrator positioned properly, the tone should be clearly audible to a person with normal hearing – less than 25 dB. When a bone vibrator fails this test, the calibration should be verified.

### Masking Level Check

Select the Tone test type. Ensure the stimulus is narrow band noise on Channel 2. Activate the Channel 2 Interrupt button and listen for a smooth, even hiss.

### Talk Forward Check

Speech should be clearly audible (in the earphones) when spoken in a normal tone with the Talk Forward dB HL control set at 45 dB HL.

## CALIBRATION

To guarantee that the AudioStar Pro works properly, the instrument should be checked and calibrated at least once per year. Have an authorized service technician perform electrical safety checks on the unit to maintain continued compliance to IEC and UL 60601-1. The use of non-calibrated devices can lead to incorrect test results.

The service and calibration must be performed by an authorized service center. If these checks are not done, Medical Device Regulation (EU) 2017/745 and other regulations may be violated, and warranties may be void.

## CLEANING THE SYSTEM

Non-disposable parts of the system, including the transducers and insert phone tubes which are in direct contact with the patient may need to be disinfected between patients.

Turn **OFF** the system and disconnect power before cleaning the instrument. Use a soft cloth lightly dampened with cleaning solution to clean all exposed surfaces. DO NOT allow liquid to contact the metal parts inside the transducers (e.g., earphones / headphones). Do not permit solutions or disinfecting agents to seep into the electronic portions of the system. Take special care around controls, connectors, and panel edges. Remove any dust from the exterior of the system with a soft brush or cloth. Use a brush to dislodge any dirt on or around the connectors and panel edges. Remove stubborn dirt with a soft cloth slightly dampened with mild detergent and water. Wipe surfaces dry afterward. Do not use instrument or transducers until they are completely dry.

Recommendations for cleaning and disinfection of the AudioStar Pro presented in this document are not intended to replace or contradict policies in effect or procedures required for infection control at the facility.

### Cleaning and Disinfecting Agents

According to the recommendations from the CDC, audiometric equipment is non-critical medical equipment and typically requires cleansing followed by low to intermediate level disinfecting, depending on the nature of the contamination. Cleaning should be done with a mild soapy detergent (such as dishwashing liquid) and a damp cloth, or an Endozime Sponge followed by an application of EPA-registered hospital disinfectant. Do not use any abrasive cleaners.

Use of a non-alcohol-based disinfectant is recommended for larger areas and headphones. Non-alcohol-based products that contain the active ingredient referred to as quaternary ammonia compound or hydrogen peroxide-based cleaner such as Oxivir Disinfectant Wipes are recommended to clean the ear cushions, headset, and to wipe down the machine. The quaternary ammonia compound and hydrogen peroxide are specifically designed to disinfect rubber, plastic, silicone, and acrylic products which are commonly used in hearing evaluation instruments.

### Single-use Accessories

Disposable or single use supplies such as eartips for the insert earphones are intended to be discarded immediately after use. Use of disposable items on multiple patients can lead to cross-contamination and infection.

## STATUS AND ERROR MESSAGES

### ERROR MESSAGES

#### Please try another selection

Indicates an incorrect selection. This could include actions such as incompatible transducers, incompatible routing, or no calibration data stored for the selected transducers.

#### No test data stored

Indicates that there is no test data available to be erased, printed or transferred.

#### Printer communication error

Occurs during printing if communication problems with the printer are encountered.

#### Record test result in comments

Test results of the ABLB and Tone Decay are not recorded directly on a report. This message indicates that the results should be documented in the comments.

#### The session comments have been updated with the results of the SDT test

This message indicates that the stored speech detection threshold results will appear in the comments section and will be printed directly or transferred electronically.

#### Not supported in speech

The selected action is not supported in the speech test type.

### Error

If there are general system errors, a dialog box with "Error" in the title will be shown with the given error.

### Error

*The startup configuration for this test type is not fully calibrated; a search for a different configuration that is calibrated has found the currently displayed configuration*

This message indicates that the selected transducers have not been calibrated.

## SUPPLIES AND ACCESSORIES

To order consumables, additional accessories and to replace detachable parts that have been damaged, please contact GSI or your GSI distributor for current prices and delivery charges. Before use check the packaging for any shelf-life expirations.

Part Number	Description
8500090	Foam Eartips for Insert Phones, 13 mm (Adult), 50 Each
8106226	Foam Eartips for Insert Phones, 13 mm (Adult), 100 Each
8500110	Foam Eartips for Insert Phones, 10 mm (Pediatric), 50 Each
8106227	Foam Eartips for Insert Phones, 10 mm (Pediatric), 100 Each
8500100	Foam Eartips for Insert Phones, Jumbo, 24 Each
8500140	Connectors (nubbins) for Insert Phones to Silicone Tubes, 10 Each
8501072	Sound tubes for Insert Phones, 1 red, 1 blue
8109155	AMBAND – Calibrated Bone Conduction Headbands (10)
8522201	Earpads for monitor headset with boom mic (2)
8504476	Microphone windscreen (Sennheiser PS 01 for PC 131 headset)
8107419	Ear cushion covers (500) 3 inch (for DD45)
8513404	Ear cushion covers (100) 4 inch (for DD450 and DD65 V2)

## APPENDIX A – TECHNICAL DATA

The AudioStar Pro is an active, diagnostic medical product. The device is classified as a class IIA device according to the Medical Device Regulation (EU) 2017/745 and a class II device according to the US FDA. Manufactured, designed, developed, and marketed under ISO 13485 certified quality systems.

<b>Standards Compliance</b>	
<b>Safety and Electromagnetic compatibility (EMC)</b>	IEC 60601-1, IEC 60601-1-2
<b>Calibration and Test Signal</b>	ANSI S3.6 (2018) Type 1 IEC 60645-1 (2017) Type 1 IEC 60645-2 (2017) ISO 389-2 ISO 389-3 ISO 389-4 ISO 389-5 ISO 389-7 ISO 389-8
<b>General Specifications</b>	
<b>ENVIRONMENTAL</b>	
<b>Transport and Handling</b> 	Transport package shall be kept away from rain and in dry conditions
<b>Temperature</b> 	Operation: + 15° C to + 40° C (+ 59° F to + 104° F) Transport: - 20° C to + 50° C (- 4° F to + 122° F) Storage: 0° C to + 50° C (+32° F to + 122° F)
<b>Humidity</b> 	Operation: Maximum relative humidity 95 %, non-condensing Transport and Storage: Maximum relative humidity 93 %, non-condensing
<b>Ambient air pressure</b>	98 kPa – 104 kPa
<b>Location</b>	Indoor Use
<b>Expected Service Life</b>	7 years

<b>POWER</b>	
<b>Power Supply</b>	100 V AC, 1A – 240 V AC, 0.6 A 50/60 Hz, 90 W
<b>PHYSICAL</b>	
<b>Dimensions</b>	51 x 37 x 33.5 cm (L x W x H) 20.1 x 14.6 x 13.2 in LCD raised
<b>Weight</b>	7.7 kg (17 lb)
<b>System</b>	
<b>Channels</b>	2, Independent
<b>Pure Tone Stimuli</b>	
<b>Frequency Range</b>	125-20,000 Hz (transducer dependent)
<b>Frequency Accuracy</b>	± 1 %
<b>Total Harmonic Distortion</b>	< 2% (earphones and insert phones) < 5.5% (B81)
<b>Level Range</b>	-20 – 120 dB HL (transducer dependent)
<b>Signal Format</b>	Steady: Tone continuously present Pulsed: Tone pulsed 200 msec ON, 200 msec OFF FM: Modulation Rate: 5 Hz, Modulation depth +/- 5% Pediatric Noise: Continuously presented or pulsed
<b>TRANSDUCERS</b>	
<b>RadioEar IP30 Insert Earphones</b>	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 120 dB HL

<b>RadioEar DD45 Headphones</b>	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 120 dB HL
<b>RadioEar DD450 Headphones</b>	Frequency Range: 125 Hz – 20,000 Hz Level: -20 dB HL to 100 dB HL
<b>RadioEar B81 Bone Conductor</b>	Frequency Range: 250 Hz – 8000 Hz Level: -10 dB HL to 90 dB HL (mastoid) Level: -10 dB HL to 80 dB HL (forehead)
<b>Sound Field</b>	Frequency Range: 125 Hz – 8000 Hz Level: -10 dB HL to 90 dB HL (basic speakers) Level: -10 dB HL to 96 dB HL (high performance speakers) Level: -10 dB HL to 102 dB HL (high performance speakers and external booster amplifier)
<b>Note:</b> The maximum HL values are applicable to the middle frequencies only	
<b>Speech</b>	
<b>Microphone</b>	For live voice testing and communications
<b>INT/EXT A &amp; INT/EXT B</b>	Can be utilized for internal speech files or recorded speech material from an external digital device
<b>Level Range</b>	Air Conduction: -10 dB HL to 100 dB HL (DD45 Linear), Inserts (IP30): -10 dB HL to 95 dB HL Bone Conduction (B81): -10 dB HL to 60 dB HL (mastoid), -10 dB HL to 50 dB HL (forehead) Sound Field: -10 dB HL to 90 dB HL (basic speakers)
<b>Masking</b>	
<b>Narrow Band Noise (effective masking)</b>	Maximum dB HL is 15 dB below tone  Frequency Range: 3 dB down bandwidth, 1/3 octave minimum; ½ octave maximum of test signal

<b>Speech Noise</b>	Air Conduction: -10 dB HL to 95 dB HL (DD45) Bone Conduction: -10 dB HL to 50 dB HL (mastoid B81), -10 dB HL to 40 dB HL (forehead B81) Sound Field: -10 dB HL to 85 dB HL  Frequency Range: Equal energy per frequency 100 to 1,000 Hz with a 12 dB/octave roll-off from 1,000 to 6,000 Hz
<b>White Noise</b>	Air Conduction: -10 dB HL to 95 dB HL (DD45) Bone Conduction: -10 dB HL to 60 dB HL (mastoid B81), -10 dB HL to 50 dB HL (forehead B81) Sound Field -10 dB HL to 80 dB HL  Frequency Range: 125 Hz to 12,000 Hz with constant bandwidth
<b>Input Specification</b>	
<b>External CD</b>	7 mV RMS at max. gain for 0 dB reading; input impedance 47 kOhm
<b>Talk Forward</b>	100 uV RMS at max. gain for 0dB reading; input impedance 3.2 kOhm
<b>Output Specification</b>	
<b>Headphones</b>	7 V RMS at 10 Ohm load; 60-20,000 Hz -3 dB
<b>Insert Phones</b>	7 V RMS at 10 Ohm load; 60-20,000 Hz -3 dB
<b>Bone</b>	7 V RMS at 10 Ohm load; 60-20,000 Hz -3 dB
<b>Sound field</b>	7 V RMS at min. 2 kOhm load; 60-20,000 kHz -3 dB
<b>Monitor</b>	2x 3 V RMS at 32 Ohm/1.5 V RMS at 8 Ohm load; 60-20,000Hz -3dB
<b>Free Field Specification</b>	

<b>Power Amplifier and Speakers</b>	With an input of 7 Vrms-Amplifier and loudspeakers must be able to create a sound pressure level of 100 dB at 1 meter and meet the following requirements
<b>Frequency Response</b>	125-250Hz, (+0/-10 dB) 250-4000Hz ( $\pm$ 3 dB) 4000-6300Hz ( $\pm$ 5 dB)
<b>Total Harmonic Distortion</b>	80 dB SPL <3% 100 dB SPL <10%
<b>Communication and Monitoring</b>	
<b>Talk Forward</b>	Permits the tester to speak through the test microphone into the selected transducer at approximately the level set by the front panel controls
<b>Talk Back</b>	Allows the tester to listen to comments from the patient in the testing booth
<b>Monitor</b>	The monitor headset or monitor speaker built into the instrument housing can be used by the tester to listen to Channel 1, Channel 2, Aux intercom, and/or Talk Back signals
<b>Aux Intercom</b>	The built-in Auxiliary Intercom and Assistant headset allows the tester to speak directly to an Assistant without the patient hearing the conversation and allows the assistant to hear what is being presented to the patient
<b>Special Tests</b>	
<b>QuickSIN</b>	Six (6) sentences with five (5) key words per sentence are presented in four-talker babble noise. The sentences are presented at pre-recorded signal-to-noise ratios. The SNR's used are 25, 20, 15, 10, 5, and 0.
<b>BKB-SIN</b>	18 List Pairs. The sentences are presented at prerecorded signal-to-noise ratios that decrease


	in 3-dB steps. Each list in the pair is individually scored, and the results of the two lists are averages to obtain the List Pair score. Results are compared to normative data to obtain the SNR Loss
<b>TEN</b>	TEN masking noise will be presented to the test ear. Pure tone stimuli between 500 and 4000 Hz can be used at 1, 2, or 5 dB increments to obtain TEN thresholds
<b>ALT (ABLB)</b>	Tone alternating between Channel 1 and Channel 2: Channel 1 is 400 msec ON, 400 msec OFF followed by Channel 2, 400 msec ON, 400 msec OFF
<b>SISI</b>	A level increment is added to a tone in the selected channel for 200 msec, every 5 seconds. The HL increments are in 1, 2 or 5 dB
<b>Others (User defined)</b>	MLB Lombard test Pure Tone Stenger Speech Stenger SAL

## APPENDIX B – SPECIFICATION OF INPUT/OUTPUT CONNECTIONS


### MAINS CABLE

Pin	Description	
1	Live	 <p>IEC C13</p>
2	Neutral	
3	Earth	


### HEADPHONES/INSERTS (2 CONNECTORS, R, L)

Pin	Description	
1	Ground	 <p>6.3mm Mono</p>
2	Signal	



### BONE OSCILLATOR

Pin	Description	
1	Ground	 <p>6.3mm Mono</p>
2	Signal	


### GOOSENECK MIC

Pin	Description	
1	Ground	 <p>6.3mm Stereo</p>
2	DC Bias	
3	Signal	


### PATIENT RESPONSE X2

Pin	Description	
1		 <p>6.3mm Mono</p>
2	Contact Closure 	


### TALK FORWARD, TALK BACK MIC

Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	DC Bias	
3	Signal	

### MONITOR HEADSET


Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	Signal 2	
3	Signal 1	

### AUX HEADSET


Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	Right	
3	Left	


### EXTERNAL A+B

Pin	Description	
1	Ground	 <p>3.5mm Stereo</p>
2	External B (CD 2)	
3	External A (CD 1)	

### VRA


Pin	Description	
1	Ground	 <p>3.5mm TRRS</p>
2	Right	
3	Center	
4	Left	

### MONITOR


Pin	Description	
1	Data2+	 <p>HDMI Type A</p>
2	Ground	
3	Data2-	
4	Data1+	
5	Ground	
6	Data1-	
7	Data0+	
8	Ground	
9	Data0-	
10	Clock+	

11	Ground
12	Clock-
13	CEC
14	Reserved
15	SCL
16	SDA
17	Ground
18	+5V DC
19	Hot Plug Detect

**SPEAKERS (FF1, FF2, FF3, FF4)**


Pin	Description	
1	Ground	 <p>RCA</p>
2	Signal	

**USB A (HOST) x4**

Pin	Description	
1	+ 5 V DC	
2	Data	
3	Data	
4	GND	

### USB B (DEVICE)

Pin	Description
1	+ 5 V DC
2	Data
3	Data
4	GND

### SPEAKER 1 &2 (5 POLE DIN)

Pin	Description
1	OUT FF+
2	Ground
3	OUT-FF2+
4	OUT-FF1-
5	OUT-FF2-

### LAN ETHERNET

Not Supported

## APPENDIX C - CALIBRATION REFERENCE & MAXIMUM LEVELS

The AudioStar Pro is supplied from the factory calibrated for the transducers that were purchased with it. The exception is the speakers, as those must be calibrated in the environment where they will be used. The calibration data supplied from the factory is only valid for GSI supplied transducers and cannot be applied to non-GSI supplied transducers.

It is recommended that calibration of the instrument and transducers be performed annually by authorized GSI Representatives using appropriate calibration instrumentation. If periodic checks are also desired, the tables in this section provide the SPL values per frequency for each transducer. If the measured values are not within  $\pm 5$  dB at 125, 6,000, and 8,000 Hz in the earphones, the GSI AudioStar Pro should be scheduled for immediate maintenance.

It is not possible to select a dB HL value outside the limits for a particular transducer/frequency combination. An attempt to change or select a hearing level control that is outside of the limit will cause the dB HL display to flash momentarily and then the test channel value will be replaced with NR (No Response). If an audiogram is displayed and the limits for a frequency/transducer are reached, the symbol for no response is displayed in the audiogram.

It is not possible to select a test frequency that is invalid for a particular transducer.

The hearing levels listed in the Max HL tables are maximum levels. These levels are achievable only if ANSI, ISO or GSI reference threshold levels, and not customized calibration values, are used. At no time will the hearing level limit exceed 120 dB HL.

## EARPHONES - PURE TONE RETSPL

Transducer	DD45	DD45	TDH50	TDH50	DD450
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
<b>125 Hz</b>	47.5	45	47.5	45	30.5
<b>160 Hz</b>	40.5	38.5	40	38.5	26
<b>200 Hz</b>	33.5	32.5	33.5	32.5	22
<b>250 Hz</b>	27	27	26.5	27	18
<b>315 Hz</b>	22.5	22	22	22	15.5
<b>400 Hz</b>	17.5	17	17.5	17	13.5
<b>500 Hz</b>	13	13.5	13.5	13.5	11
<b>630 Hz</b>	9	10.5	10.5	10.5	8
<b>750 Hz</b>	6.5	9	8.5	9	6
<b>800 Hz</b>	6.5	8.5	8.5	8.5	6
<b>1000 Hz</b>	6	7.5	7.5	7.5	5.5
<b>1250 Hz</b>	7	7.5	7.5	7.5	6
<b>1500 Hz</b>	8	7.5	7.5	7.5	5.5
<b>1600 Hz</b>	8	8	8.5	8	5.5
<b>2000 Hz</b>	8	9	11	9	4.5
<b>2500 Hz</b>	8	10.5	10	10.5	3
<b>3000 Hz</b>	8	11.5	9.5	11.5	2.5
<b>3150 Hz</b>	8	11.5	9.5	11.5	4
<b>4000 Hz</b>	9	12	10.5	12	9.5
<b>5000 Hz</b>	13	11	12	11	14
<b>6000 Hz</b>	20.5	16	13.5	16	17
<b>6300 Hz</b>	19	21	13.5	21	17.5
<b>8000 Hz</b>	12	15.5	13	15.5	17.5
<b>9000 Hz</b>					19
<b>10000 Hz</b>					22
<b>11200 Hz</b>					23
<b>12000 Hz</b>			17.5	11	0
<b>12500 Hz</b>					27.5
<b>14000 Hz</b>					35
<b>16000 Hz</b>					56
<b>18000 Hz</b>					83
<b>20000 Hz</b>					105

- DD45 6ccm uses IEC60318-3 or NBS 9A coupler and RETSPL comes from PTB – DTU report 2009-2010. Force 4.5N  $\pm$ 0.5N
- DD45 Artificial ear uses IEC60318-1 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-1 1998. Force 4.5N  $\pm$ 0.5N
- TDH50 6ccm uses IEC60318-3 or NBS 9A coupler and RETSPL comes from ANSI S3.6 2018. Force 4.5N  $\pm$ 0.5N
- TDH50 Artificial ear uses IEC60318-1 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-1 1998 Force 4.5N  $\pm$ 0.5N
- HDA200 Artificial ear uses IEC60318-1 coupler with type 1 adaptor and RETSPL comes from ANSI S3.6 2018 and ISO 389-8 2004. Force 9N  $\pm$ 0.5N

## EARPHONES - ANSI SPEECH RETSPL

Transducer	DD45	DD45	TDH50	TDH50	DD450
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
<b>Speech</b>	18.5	20	20	20	19
<b>Speech Equ.FF.</b>	18.5	19.5	17	18	18.5
<b>Speech Non-linear</b>	6	7.5	7.5	7.5	5.5
<b>Speech noise</b>	18.5	20	20	20	19
<b>Speech noise Equ.FF.</b>	18.5	19.5	17	18	18.5
<b>Speech noise Non-linear</b>	6	7.5	7.5	7.5	5.5
<b>White noise in speech</b>	21	22.5	22.5	22.5	21.5

- DD45 (G<sub>F</sub>-G<sub>C</sub>) PTB-DTU report 2009-2010.
- TDH50 (G<sub>F</sub>-G<sub>C</sub>) ANSI S3.6 2018.
- HDA200 (G<sub>F</sub>-G<sub>C</sub>) ANSI S3.6 2018 and ISO 389-8 2004.
- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL - (G<sub>F</sub>-G<sub>C</sub>) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

## EARPHONES - IEC SPEECH RETSPL

<b>Transducer</b>	<b>DD45</b>	<b>DD45</b>	<b>TDH50</b>	<b>TDH50</b>	<b>DD450</b>
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	RETSPL	RETSPL	RETSPL	RETSPL	RETSPL
<b>Speech</b>	20	20	20	20	20
<b>Speech Equ.FF.</b>	3.5	4.5	2	3	3.5
<b>Speech Non-linear</b>	6	7.5	7.5	7.5	5.5
<b>Speech noise</b>	20	20	20	20	20
<b>Speech noise Equ.FF.</b>	3.5	4.5	2	3	3.5
<b>Speech noise Non-linear</b>	6	7.5	7.5	7.5	5.5
<b>White noise in speech</b>	22.5	22.5	22.5	22.5	22.5

- DD45 (G<sub>F</sub>-G<sub>C</sub>) PTB-DTU report 2009-2010.
- TDH50 (G<sub>F</sub>-G<sub>C</sub>) ANSI S3.6 2018.
- HDA200 (G<sub>F</sub>-G<sub>C</sub>) ANSI S3.6 2018 and ISO 389-8 2004.
- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level (G<sub>F</sub>-G<sub>C</sub>) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 - B81 IEC60645-2 2017 (no weighting)

## EARPHONES - PURE TONE MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
<b>Tone 125 Hz</b>	90	90	85	85	100
<b>Tone 160 Hz</b>	95	95	90	90	105
<b>Tone 200 Hz</b>	100	100	95	95	105
<b>Tone 250 Hz</b>	110	110	105	105	110
<b>Tone 315 Hz</b>	115	115	110	110	115
<b>Tone 400 Hz</b>	120	120	115	115	115
<b>Tone 500 Hz</b>	120	120	120	120	115
<b>Tone 630 Hz</b>	120	120	120	120	120
<b>Tone 750 Hz</b>	120	120	120	120	120
<b>Tone 800 Hz</b>	120	120	120	120	120
<b>Tone 1000 Hz</b>	120	120	120	120	120
<b>Tone 1250 Hz</b>	120	120	120	120	110
<b>Tone 1500 Hz</b>	120	120	120	120	115
<b>Tone 1600 Hz</b>	120	120	120	120	115
<b>Tone 2000 Hz</b>	120	120	120	120	115
<b>Tone 2500 Hz</b>	120	120	120	120	115
<b>Tone 3000 Hz</b>	120	120	120	120	115
<b>Tone 3150 Hz</b>	120	120	120	120	115
<b>Tone 4000 Hz</b>	120	120	120	120	115
<b>Tone 5000 Hz</b>	120	120	115	115	105
<b>Tone 6000 Hz</b>	115	115	115	110	105
<b>Tone 6300 Hz</b>	115	110	110	105	105
<b>Tone 8000 Hz</b>	110	110	100	100	105
<b>Tone 9000 Hz</b>					100
<b>Tone 10000 Hz</b>					100
<b>Tone 11200 Hz</b>					95
<b>Tone 12000 Hz</b>			90	90	
<b>Tone 12500 Hz</b>					90
<b>Tone 14000 Hz</b>					80
<b>Tone 16000 Hz</b>					60
<b>Tone 18000 Hz</b>					30
<b>Tone 20000 Hz</b>					15

## EARPHONES - NB NOISE EFFECTIVE MASKING LEVEL

<b>Transducer</b>	<b>DD45</b>	<b>DD45</b>	<b>TDH50</b>	<b>TDH50</b>	<b>DD450</b>
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	EM	EM	EM	EM	EM
<b>NB 125 Hz</b>	51.5	49	51.5	49	34.5
<b>NB 160 Hz</b>	44.5	42.5	44	42.5	30
<b>NB 200 Hz</b>	37.5	36.5	37.5	36.5	26
<b>NB 250 Hz</b>	31	31	30.5	31	22
<b>NB 315 Hz</b>	26.5	26	26	26	19.5
<b>NB 400 Hz</b>	21.5	21	21.5	21	17.5
<b>NB 500 Hz</b>	17	17.5	17.5	17.5	15
<b>NB 630 Hz</b>	14	15.5	15.5	15.5	13
<b>NB 750 Hz</b>	11.5	14	13.5	14	11
<b>NB 800 Hz</b>	11.5	13.5	13.5	13.5	11
<b>NB 1000 Hz</b>	12	13.5	13.5	13.5	11.5
<b>NB 1250 Hz</b>	13	13.5	13.5	13.5	12
<b>NB 1500 Hz</b>	14	13.5	13.5	13.5	11.5
<b>NB 1600 Hz</b>	14	14	14.5	14	11.5
<b>NB 2000 Hz</b>	14	15	17	15	10.5
<b>NB 2500 Hz</b>	14	16.5	16	16.5	9
<b>NB 3000 Hz</b>	14	17.5	15.5	17.5	8.5
<b>NB 3150 Hz</b>	14	17.5	15.5	17.5	10
<b>NB 4000 Hz</b>	14	17	15.5	17	14.5
<b>NB 5000 Hz</b>	18	16	17	16	19
<b>NB 6000 Hz</b>	25.5	21	18.5	21	22
<b>NB 6300 Hz</b>	24	26	18.5	26	22.5
<b>NB 8000 Hz</b>	17	20.5	18	20.5	22.5
<b>NB 9000 Hz</b>					24
<b>NB 10000 Hz</b>					27
<b>NB 11200 Hz</b>					28
<b>NB 12000 Hz</b>			22.5	16	
<b>NB 12500 Hz</b>					32.5
<b>NB 14000 Hz</b>					40
<b>NB 16000 Hz</b>					61
<b>NB 18000 Hz</b>					88
<b>NB 20000 Hz</b>					110
<b>White noise</b>	0	0	0	0	0
<b>TEN noise</b>	25	25	24.5	24.5	31.5

## EARPHONES - NB NOISE MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
Impedance	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
Coupler	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
<b>NB 125 Hz</b>	75	75	65	65	75
<b>NB 160 Hz</b>	80	80	70	70	80
<b>NB 200 Hz</b>	90	90	80	80	80
<b>NB 250 Hz</b>	95	95	85	85	85
<b>NB 315 Hz</b>	100	100	90	90	90
<b>NB 400 Hz</b>	105	105	95	95	95
<b>NB 500 Hz</b>	110	110	100	100	95
<b>NB 630 Hz</b>	110	110	100	100	95
<b>NB 750 Hz</b>	110	110	105	105	100
<b>NB 800 Hz</b>	110	110	105	105	100
<b>NB 1000 Hz</b>	110	110	105	105	100
<b>NB 1250 Hz</b>	110	110	105	105	95
<b>NB 1500 Hz</b>	110	110	105	105	100
<b>NB 1600 Hz</b>	110	110	105	105	100
<b>NB 2000 Hz</b>	110	110	100	100	100
<b>NB 2500 Hz</b>	110	110	100	100	100
<b>NB 3000 Hz</b>	110	110	100	100	100
<b>NB 3150 Hz</b>	110	110	100	100	100
<b>NB 4000 Hz</b>	110	110	100	100	100
<b>NB 5000 Hz</b>	110	110	100	100	95
<b>NB 6000 Hz</b>	105	105	95	95	90
<b>NB 6300 Hz</b>	105	100	95	90	90
<b>NB 8000 Hz</b>	100	100	90	85	90
<b>NB 9000 Hz</b>					85
<b>NB 10000 Hz</b>					85
<b>NB 11200 Hz</b>					80
<b>NB 12000 Hz</b>			75	75	
<b>NB 12500 Hz</b>					75
<b>NB 14000 Hz</b>					70
<b>NB 16000 Hz</b>					50
<b>NB 18000 Hz</b>					20
<b>NB 20000 Hz</b>					0
<b>White noise</b>	120	120	120	120	115
<b>TEN noise</b>	110	110	100	100	121

## EARPHONES - ANSI SPEECH MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
<b>Speech</b>	110	105	100	100	90
<b>Speech Equ.FF.</b>	100	100	90	90	85
<b>Speech Non-linear</b>	120	120	115	115	110
<b>Speech noise</b>	100	100	95	90	85
<b>Speech noise Equ.FF.</b>	100	95	85	85	80
<b>Speech noise Non-linear</b>	115	115	110	110	105
<b>White noise in speech</b>	95	95	95	95	90

## EARPHONES - IEC SPEECH MAX HL

Transducer	DD45	DD45	TDH50	TDH50	DD450
<b>Impedance</b>	10 $\Omega$	10 $\Omega$	60 $\Omega$	60 $\Omega$	23 $\Omega$
<b>Coupler</b>	318-3	318-1	318-3	318-1	318-1
	Max HL	Max HL	Max HL	Max HL	Max HL
<b>Speech</b>	110	105	100	100	90
<b>Speech Equ.FF.</b>	115	115	105	105	100
<b>Speech Non-linear</b>	120	120	115	115	110
<b>Speech noise</b>	100	100	95	90	85
<b>Speech noise Equ.FF.</b>	115	110	100	100	95
<b>Speech noise Non-linear</b>	115	115	110	110	105
<b>White noise in speech</b>	95	95	95	95	90

INSERT EARPHONES - PURE TONE RETSPL

Transducer	IP30	IP30
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	RETSPL	RETSPL
<b>125 Hz</b>	26	28
<b>160 Hz</b>	22	24.5
<b>200 Hz</b>	18	21.5
<b>250 Hz</b>	14	17.5
<b>315 Hz</b>	12	15.5
<b>400 Hz</b>	9	13
<b>500 Hz</b>	5.5	9.5
<b>630 Hz</b>	4	7.5
<b>750 Hz</b>	2	6
<b>800 Hz</b>	1.5	5.5
<b>1000 Hz</b>	0	5.5
<b>1250 Hz</b>	2	8.5
<b>1500 Hz</b>	2	9.5
<b>1600 Hz</b>	2	9.5
<b>2000 Hz</b>	3	11.5
<b>2500 Hz</b>	5	13.5
<b>3000 Hz</b>	3.5	13
<b>3150 Hz</b>	4	13
<b>4000 Hz</b>	5.5	15
<b>5000 Hz</b>	5	18.5
<b>6000 Hz</b>	2	16
<b>6300 Hz</b>	2	16
<b>8000 Hz</b>	0	15.5
<b>9000 Hz</b>		
<b>10000 Hz</b>		
<b>11200 Hz</b>		
<b>12000 Hz</b>		
<b>12500 Hz</b>		
<b>14000 Hz</b>		
<b>16000 Hz</b>		
<b>18000 Hz</b>		
<b>20000 Hz</b>		

- IP30 2ccm uses ANSI S3.7-1995 IEC60318-5 coupler (HA-2 with 5mm rigid Tube) and RETSPL comes from ANSI S3.6 2018 and ISO 389-2 1994
- IP30 Ear simulator uses ANSI S3.25- IEC60318-4 coupler and RETSPL comes from ANSI S3.6 2018 and ISO 389-2 1994

## INSERT EARPHONES - ANSI SPEECH RETSPL

Transducer	IP30	IP30
<b>Impedance</b>	10 Ω / 50 Ω	10 Ω / 50 Ω
<b>Coupler</b>	2ccm	711
	RETSPL	RETSPL
<b>Speech</b>		
<b>Speech Equ.FF.</b>		
<b>Speech Non-linear</b>	12.5	18
<b>Speech noise</b>		
<b>Speech noise Equ.FF.</b>		
<b>Speech noise Non-linear</b>	12.5	18
<b>White noise in speech</b>	15	20.5

- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL – (GF-GC) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

## INSERT EARPHONES - IEC SPEECH RETSPL

Transducer	IP30	IP30
<b>Impedance</b>	10 Ω / 50 Ω	10 Ω / 50 Ω
<b>Coupler</b>	2ccm	711
	RETSPL	RETSPL
<b>Speech</b>		
<b>Speech Equ.FF.</b>		
<b>Speech Non-linear</b>	20	20
<b>Speech noise</b>		
<b>Speech noise Equ.FF.</b>		
<b>Speech noise Non-linear</b>	20	20
<b>White noise in speech</b>	22.5	22.5

- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level ( $G_F-G_C$ ) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 -B81 IEC60645-2 2017 (no weighting)

INSERT EARPHONES - PURE TONE MAX HL

<b>Transducer</b>	<b>IP30</b>	<b>IP30</b>
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	Max HL	Max HL
<b>Tone 125 Hz</b>	90	90
<b>Tone 160 Hz</b>	95	95
<b>Tone 200 Hz</b>	100	100
<b>Tone 250 Hz</b>	105	105
<b>Tone 315 Hz</b>	105	105
<b>Tone 400 Hz</b>	110	110
<b>Tone 500 Hz</b>	110	110
<b>Tone 630 Hz</b>	115	115
<b>Tone 750 Hz</b>	115	115
<b>Tone 800 Hz</b>	115	115
<b>Tone 1000 Hz</b>	120	120
<b>Tone 1250 Hz</b>	120	120
<b>Tone 1500 Hz</b>	120	120
<b>Tone 1600 Hz</b>	120	120
<b>Tone 2000 Hz</b>	120	120
<b>Tone 2500 Hz</b>	120	120
<b>Tone 3000 Hz</b>	120	120
<b>Tone 3150 Hz</b>	120	120
<b>Tone 4000 Hz</b>	115	115
<b>Tone 5000 Hz</b>	105	105
<b>Tone 6000 Hz</b>	100	100
<b>Tone 6300 Hz</b>	100	100
<b>Tone 8000 Hz</b>	95	95
<b>Tone 9000 Hz</b>		
<b>Tone 10000 Hz</b>		
<b>Tone 11200 Hz</b>		
<b>Tone 12000 Hz</b>		
<b>Tone 12500 Hz</b>		
<b>Tone 14000 Hz</b>		
<b>Tone 16000 Hz</b>		
<b>Tone 18000 Hz</b>		
<b>Tone 20000 Hz</b>		

INSERT EARPHONES - NB NOISE EFFECTIVE MASKING LEVEL

Transducer	IP30	IP30
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	EM	EM
<b>NB 125 Hz</b>	30	32
<b>NB 160 Hz</b>	26	28.5
<b>NB 200 Hz</b>	22	25.5
<b>NB 250 Hz</b>	18	21.5
<b>NB 315 Hz</b>	16	19.5
<b>NB 400 Hz</b>	13	17
<b>NB 500 Hz</b>	9.5	13.5
<b>NB 630 Hz</b>	9	12.5
<b>NB 750 Hz</b>	7	11
<b>NB 800 Hz</b>	6.5	10.5
<b>NB 1000 Hz</b>	6	11.5
<b>NB 1250 Hz</b>	8	14.5
<b>NB 1500 Hz</b>	8	15.5
<b>NB 1600 Hz</b>	8	15.5
<b>NB 2000 Hz</b>	9	17.5
<b>NB 2500 Hz</b>	11	19.5
<b>NB 3000 Hz</b>	9.5	19
<b>NB 3150 Hz</b>	10	19
<b>NB 4000 Hz</b>	10.5	20
<b>NB 5000 Hz</b>	10	23.5
<b>NB 6000 Hz</b>	7	21
<b>NB 6300 Hz</b>	7	21
<b>NB 8000 Hz</b>	5	20.5
<b>NB 9000 Hz</b>		
<b>NB 10000 Hz</b>		
<b>NB 11200 Hz</b>		
<b>NB 12000 Hz</b>		
<b>NB 12500 Hz</b>		
<b>NB 14000 Hz</b>		
<b>NB 16000 Hz</b>		
<b>NB 18000 Hz</b>		
<b>NB 20000 Hz</b>		
<b>White noise</b>	0	0
<b>TEN noise</b>	16	25

INSERT EARPHONES - NB NOISE MAX HL

Transducer	IP30	IP30
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	Max HL	Max HL
<b>NB 125 Hz</b>	90	90
<b>NB 160 Hz</b>	95	95
<b>NB 200 Hz</b>	100	100
<b>NB 250 Hz</b>	105	105
<b>NB 315 Hz</b>	105	105
<b>NB 400 Hz</b>	105	105
<b>NB 500 Hz</b>	110	110
<b>NB 630 Hz</b>	110	110
<b>NB 750 Hz</b>	110	110
<b>NB 800 Hz</b>	110	110
<b>NB 1000 Hz</b>	110	110
<b>NB 1250 Hz</b>	110	110
<b>NB 1500 Hz</b>	110	110
<b>NB 1600 Hz</b>	110	110
<b>NB 2000 Hz</b>	110	110
<b>NB 2500 Hz</b>	110	110
<b>NB 3000 Hz</b>	110	110
<b>NB 3150 Hz</b>	110	110
<b>NB 4000 Hz</b>	110	110
<b>NB 5000 Hz</b>	105	105
<b>NB 6000 Hz</b>	100	100
<b>NB 6300 Hz</b>	100	100
<b>NB 8000 Hz</b>	95	90
<b>NB 9000 Hz</b>		
<b>NB 10000 Hz</b>		
<b>NB 11200 Hz</b>		
<b>NB 12000 Hz</b>		
<b>NB 12500 Hz</b>		
<b>NB 14000 Hz</b>		
<b>NB 16000 Hz</b>		
<b>NB 18000 Hz</b>		
<b>NB 20000 Hz</b>		
<b>White noise</b>	110	110
<b>TEN noise</b>	100	100

INSERT EARPHONES - ANSI SPEECH MAX HL

<b>Transducer</b>	<b>IP30</b>	<b>IP30</b>
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	Max HL	Max HL
<b>Speech (Non-linear)</b>	110	105
<b>Speech noise Non-linear</b>	100	100
<b>White noise in speech</b>	95	85

INSERT EARPHONES - IEC SPEECH MAX HL

<b>Transducer</b>	<b>IP30</b>	<b>IP30</b>
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	2ccm	711
	Max HL	Max HL
<b>Speech (Non-linear)</b>	100	105
<b>Speech noise Non-linear</b>	90	100
<b>White noise in speech</b>	85	85

BONE VIBRATOR - PURE TONE RETFL

<b>Transducer</b>	<b>B81</b>	<b>B81</b>
<b>Impedance</b>	10 $\Omega$	10 $\Omega$
<b>Coupler</b>	Mastoid	Forehead
	RETFL	RETFL
<b>125 Hz</b>		
<b>160 Hz</b>		
<b>200 Hz</b>		
<b>250 Hz</b>	67	79
<b>315 Hz</b>	64	76.5
<b>400 Hz</b>	61	74.5
<b>500 Hz</b>	58	72
<b>630 Hz</b>	52.5	66
<b>750 Hz</b>	48.5	61.5
<b>800 Hz</b>	47	59
<b>1000 Hz</b>	42.5	51
<b>1250 Hz</b>	39	49
<b>1500 Hz</b>	36.5	47.5
<b>1600 Hz</b>	35.5	46.5
<b>2000 Hz</b>	31	42.5
<b>2500 Hz</b>	29.5	41.5
<b>3000 Hz</b>	30	42
<b>3150 Hz</b>	31	42.5
<b>4000 Hz</b>	35.5	43.5
<b>5000 Hz</b>	40	51
<b>6000 Hz</b>	40	51
<b>6300 Hz</b>	40	50
<b>8000 Hz</b>	40	50
<b>9000 Hz</b>		
<b>10000 Hz</b>		
<b>11200 Hz</b>		
<b>12000 Hz</b>		
<b>12500 Hz</b>		
<b>14000 Hz</b>		
<b>16000 Hz</b>		
<b>18000 Hz</b>		
<b>20000 Hz</b>		

- B81 uses ANSI S3.13 or IEC60318-6 2007 mechanical coupler and RETFL come from ANSI S3.6 2018 and ISO 389-3 1994. Force 5.4N  $\pm$ 0.5N

## BONE VIBRATOR - ANSI SPEECH RETSPL

Transducer	B81	B81
Impedance	10 $\Omega$	10 $\Omega$
Coupler	Mastoid	Forehead
	RETFL	RETFL
Speech (Non-linear)	55	63.5
Speech noise Non-linear	55	63.5
White noise in speech	57.5	66

- ANSI Speech level 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (acoustical linear weighting)
- ANSI Speech Equivalent free field level 12.5 dB + 1 kHz RETSPL - (GF-GC) from ANSI S3.6 2018(acoustical equivalent sensitivity weighting)
- ANSI Speech Not linear level 1 kHz RETSPL ANSI S3.6 2018 (DD45-TDH50-HDA200) and IP30- B81 12.5 dB + 1 kHz RETSPL ANSI S3.6 2018 (no weighting)

## BONE VIBRATOR - IEC SPEECH RETSPL

Transducer	B81	B81
Impedance	10 $\Omega$	10 $\Omega$
Coupler	Mastoid	Forehead
	RETFL	RETFL
Speech (Non-linear)	55	63.5
Speech noise Non-linear	55	63.5
White noise in speech	57.5	66

- IEC Speech level IEC60645-2 2017 (acoustical linear weighting)
- IEC Speech Equivalent free field level ( $G_F-G_C$ ) from IEC60645-2 2017 (acoustical equivalent sensitivity weighting)
- IEC Speech Not linear level 1 kHz RETSPL (DD45-TDH50-HDA200) and IP30 - B81 IEC60645-2 2017 (no weighting)

BONE VIBRATOR - PURE TONE MAX HL

Transducer	B81	B81
Impedance	10 Ω	10 Ω
Coupler	Mastoid	Forehead
	Max HL	Max HL
Tone 125 Hz		
Tone 160 Hz		
Tone 200 Hz		
Tone 250 Hz	50	35
Tone 315 Hz	60	45
Tone 400 Hz	70	55
Tone 500 Hz	70	55
Tone 630 Hz	75	60
Tone 750 Hz	75	60
Tone 800 Hz	75	60
Tone 1000 Hz	85	75
Tone 1250 Hz	90	80
Tone 1500 Hz	90	80
Tone 1600 Hz	90	75
Tone 2000 Hz	90	75
Tone 2500 Hz	85	70
Tone 3000 Hz	85	70
Tone 3150 Hz	85	70
Tone 4000 Hz	85	70
Tone 5000 Hz	70	55
Tone 6000 Hz	60	50
Tone 6300 Hz	55	45
Tone 8000 Hz	50	40
Tone 9000 Hz		
Tone 10000 Hz		
Tone 11200 Hz		
Tone 12000 Hz		
Tone 12500 Hz		
Tone 14000 Hz		
Tone 16000 Hz		
Tone 18000 Hz		
Tone 20000 Hz		

BONE VIBRATOR - NB NOISE EFFECTIVE MASKING LEVEL

<b>Transducer</b>	<b>B81</b>	<b>B81</b>
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	Mastoid	Forehead
	EM	EM
<b>NB 125 Hz</b>		
<b>NB 160 Hz</b>		
<b>NB 200 Hz</b>		
<b>NB 250 Hz</b>	71	83
<b>NB 315 Hz</b>	68	80.5
<b>NB 400 Hz</b>	65	78.5
<b>NB 500 Hz</b>	62	76
<b>NB 630 Hz</b>	57.5	71
<b>NB 750 Hz</b>	53.5	66.5
<b>NB 800 Hz</b>	52	64
<b>NB 1000 Hz</b>	48.5	57
<b>NB 1250 Hz</b>	45	55
<b>NB 1500 Hz</b>	42.5	53.5
<b>NB 1600 Hz</b>	41.5	52.5
<b>NB 2000 Hz</b>	37	48.5
<b>NB 2500 Hz</b>	35.5	47.5
<b>NB 3000 Hz</b>	36	48
<b>NB 3150 Hz</b>	37	48.5
<b>NB 4000 Hz</b>	40.5	48.5
<b>NB 5000 Hz</b>	45	56
<b>NB 6000 Hz</b>	45	56
<b>NB 6300 Hz</b>	45	55
<b>NB 8000 Hz</b>	45	55
<b>NB 9000 Hz</b>		
<b>NB 10000 Hz</b>		
<b>NB 11200 Hz</b>		
<b>NB 12000 Hz</b>		
<b>NB 12500 Hz</b>		
<b>NB 14000 Hz</b>		
<b>NB 16000 Hz</b>		
<b>NB 18000 Hz</b>		
<b>NB 20000 Hz</b>		
<b>White noise</b>	42.5	51
<b>TEN noise</b>		

BONE VIBRATOR - NB NOISE MAX HL

<b>Transducer</b>	<b>B81</b>	<b>B81</b>
<b>Impedance</b>	10 Ω	10 Ω
<b>Coupler</b>	Mastoid	Forehead
	Max HL	Max HL
<b>NB 125 Hz</b>		
<b>NB 160 Hz</b>		
<b>NB 200 Hz</b>		
<b>NB 250 Hz</b>	40	25
<b>NB 315 Hz</b>	50	35
<b>NB 400 Hz</b>	60	45
<b>NB 500 Hz</b>	60	45
<b>NB 630 Hz</b>	65	50
<b>NB 750 Hz</b>	65	50
<b>NB 800 Hz</b>	65	50
<b>NB 1000 Hz</b>	70	60
<b>NB 1250 Hz</b>	75	60
<b>NB 1500 Hz</b>	75	60
<b>NB 1600 Hz</b>	75	60
<b>NB 2000 Hz</b>	70	55
<b>NB 2500 Hz</b>	65	50
<b>NB 3000 Hz</b>	65	50
<b>NB 3150 Hz</b>	65	50
<b>NB 4000 Hz</b>	60	50
<b>NB 5000 Hz</b>	55	45
<b>NB 6000 Hz</b>	50	40
<b>NB 6300 Hz</b>	45	35
<b>NB 8000 Hz</b>	40	30
<b>NB 9000 Hz</b>		
<b>NB 10000 Hz</b>		
<b>NB 11200 Hz</b>		
<b>NB 12000 Hz</b>		
<b>NB 12500 Hz</b>		
<b>NB 14000 Hz</b>		
<b>NB 16000 Hz</b>		
<b>NB 18000 Hz</b>		
<b>NB 20000 Hz</b>		
<b>White noise</b>	70	60
<b>TEN noise</b>		

## BONE VIBRATOR - ANSI SPEECH MAX HL

<b>Transducer</b>	<b>B81</b>	<b>B81</b>
<b>Impedance</b>	10 $\Omega$	10 $\Omega$
<b>Coupler</b>	Mastoid	Forehead
	Max HL	Max HL
<b>Speech</b>		
<b>Speech Equ.FF.</b>		
<b>Speech Non-linear</b>	60	50
<b>Speech noise</b>		
<b>Speech noise Equ.FF.</b>		
<b>Speech noise Non-linear</b>	50	40
<b>White noise in speech</b>	60	50

## BONE VIBRATOR - IEC SPEECH MAX HL

<b>Transducer</b>	<b>B81</b>	<b>B81</b>
<b>Impedance</b>	10 $\Omega$	10 $\Omega$
<b>Coupler</b>	Mastoid	Forehead
	Max HL	Max HL
<b>Speech</b>		
<b>Speech Equ.FF.</b>		
<b>Speech Non-linear</b>	60	50
<b>Speech noise</b>		
<b>Speech noise Equ.FF.</b>		
<b>Speech noise Non-linear</b>	50	40
<b>White noise in speech</b>	60	50

FREE FIELD SPEAKERS – ANSI RETSPL AND MAX HL

Hz	Binaural RETSPL dB					Max HL dB <sup>1</sup>		
						Basic Speakers	GSI High Performance Speakers	GSI High Performance Speakers External Amplifier
	0° <sup>4</sup>	45° <sup>4</sup>	90° <sup>4</sup>	135° <sup>3</sup>	180° <sup>3</sup>	45°	45°	45°
<b>125</b>	22	21.5	21	21.5	22	70	70	75
<b>160</b>	18	17	16.5	17	18	75	75	75
<b>200</b>	14.5	13.5	13	13.5	14.5	80	80	75
<b>250</b>	11.5	10.5	9.5	10	12	80	80	90
<b>315</b>	8.5	7	6	7	9	80	80	90
<b>400</b>	6	3.5	2.5	3.5	6.5	90	90	95
<b>500</b>	4.5	1.5	0	1.5	5	95	100	110
<b>630</b>	3	-0.5	-2	-0.5	3.5	95	100	110
<b>750</b>	2.5	-1	-2.5	-1.5	2.5	95	105	110
<b>800</b>	2	-1.5	-3	-2	1.5	95	105	110
<b>1000</b>	2.5	-1.5	-3	-2.5	0.5	95	105	110
<b>1250</b>	3.5	-0.5	-2.5	-2.5	1	95	105	110
<b>1500</b>	2.5	1	-2.5	-2	1	100	105	110
<b>1600</b>	1.5	-2	-3	-2.5	0	100	105	110
<b>2000</b>	-1.5	-4.5	-3.5	-2	0	100	105	115
<b>2500</b>	-4	-7.5	-6	-2	-1.5	100	105	115
<b>3000</b>	-6	-11	-8.5	-4	-3.5	105	105	115
<b>3150</b>	-6	-11	-8	-4	-3.5	105	105	115
<b>4000</b>	-5.5	-9.5	-5	-0.5	-0.5	105	105	115
<b>5000</b>	-1.5	-7.5	-5.5	2.5	3	100	100	105
<b>6000</b>	4.5	-3	-5	3	9.5	95	95	105
<b>6300</b>	6	-1.5	-4	2.5	11	95	95	100
<b>8000</b>	12.5	7	4.1	8.5	16.5	80	80	75
<b>9000</b>	14	9	7	10.5	18.5	85	85	85
<b>10000</b>	14	9.5	8	10	18	80	80	80
<b>11200</b>	13	9	6	8	16	80	80	80
<b>12500</b>	12.5	11	4.5	6.5	15	75	75	75
<b>Speech</b>	14.5	12.5	11	10	13	90 <sup>2</sup>	96 <sup>2</sup>	102 <sup>2</sup>

- 1) When using the Calibration App, the Custom Max procedure must be utilized.
- 2) To achieve the maximum output level the crest factor per ANSI S3.6 (2018) is reduced.
- 3) The angle 135 and 180 degrees are from E.A.G Shaw and M.M Vaillancourt publication 1985
- 4) Based on ANSI S3.6-2018 / ISO 389-7

## APPENDIX D - WARRANTY AND REPAIR

We, Grason-Stadler, warrant that this product is free from defects in material and workmanship and, when properly installed and used, will perform in accordance with applicable specifications. This warranty is extended to the original purchaser of the instrument by GSI through the distributor from whom it was purchased and covers defects in material and workmanship for a period of one year from date of delivery of the instrument to the original purchaser. If within one year after original shipment, it is found not to meet this standard; it will be repaired, or at our option, replaced at no charge except for transportation costs, when returned to an authorized Grason-Stadler facility. If field service is requested, there will be no charge for labor or material; however, there will be a charge for travel expense at the service center's current rate.

---

**NOTE:** Opening the instrument case or changes to the product not approved in writing by Grason-Stadler shall void this warranty. Grason-Stadler shall not be responsible for any indirect, special, or consequential damages, even if notice has been given in advance of the possibility of such damage. The transducers may go out of calibration due to rough handling or impact (dropping).

---

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

We advise you against attempting to rectify any faults by yourself or commissioning non-experts to do so. Equipment is not user repairable. Repairs must be performed by an authorized service representative only.

No modifications of the equipment are allowed by anyone other than a qualified GSI representative. Modification of the equipment could be hazardous. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

To ensure that your instrument works properly, the GSI AudioStar Pro should be checked and calibrated at least once per year. This check must be carried out by your dealer or authorized GSI service facility.

When returning the instrument for repairs or calibration it is essential to send the acoustic transducers with the device. Send the device to an authorized service center only. Please include a detailed description of faults. To prevent damage in transit, please use the original packing material, if possible, when returning the instrument.

## APPENDIX E – RECYCLING AND DISPOSAL

Many local laws and regulations require special procedures to recycle or dispose of electrical equipment and related waste including batteries, printed circuit boards, electronic components, wiring and other elements of electronic devices. Follow all local laws and regulations for the proper disposal of batteries and any other parts of this system.

Below is the contact address for proper return or disposal of electronic wastes relating to Grason-Stadler products in Europe and other localities. The contact information for the WEEE in Europe:

Grason-Stadler  
c/o DGS Diagnostics A/S  
Audiometer Alle 1  
5500 Middelfart  
Denmark

Within the European Union, it is illegal to dispose of electric and electronic items in unsorted municipal waste. Electric and electronic waste may contain hazardous substances and therefore must be collected separately. Such products will be marked with the crossed-out wheeled bin symbol. The cooperation of the user is important to ensure a high level of reuse and recycling of electric and electronic waste. Failing to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

## APPENDIX F – ELECTROMAGNETIC COMPATIBILITY (EMC)

Portable and mobile RF communications equipment can affect the AudioStar Pro. Install and operate the system according to the EMC information presented in this chapter.

The AudioStar Pro has been tested for EMC emissions and immunity as a stand-alone device. Do not use the system adjacent to or stacked with other electronic equipment. If adjacent or stacked use is necessary, the user should verify normal operation in the configuration.

The use of accessories, transducers, and cables other than those specified, except for servicing parts sold by GSI as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the device.

Anyone connecting additional equipment is responsible for making sure the system complies with the IEC 60601-1-2 standard.

### CAUTIONS REGARDING EMC

This instrument is suitable in hospital environments except for near active HF surgical equipment and RF shielded rooms of systems for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.

Use of this instrument adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this instrument and the other equipment should be observed to verify that they are operating normally.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the AudioStar Pro, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The AudioStar Pro is intended for use in the electromagnetic environment specified below. The customer or the user of the Instrument should assure that it is used in such an environment.

**GUIDANCE AND MANUFACTURER'S DECLARATION  
ELECTROMAGNETIC EMISSIONS**

<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B Limits	The AudioStar Pro is suitable for use in all commercial, industrial, business, and residential environments.
Harmonic emissions IEC 61000-3-2	Class A Category	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

## RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT

The AudioStar Pro is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.17\sqrt{P}$	80 MHz to 800 MHz $d = 1.17\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.23\sqrt{P}$
<b>0.01</b>	0.12	0.12	0.23
<b>0.1</b>	0.37	0.37	0.74
<b>1</b>	1.17	1.17	2.33
<b>10</b>	3.70	3.70	7.37
<b>100</b>	11.70	11.70	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

**Note 1** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2** These guidelines may not apply to all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## GUIDANCE AND MANUFACTURER'S DECLARATION

### ELECTROMAGNETIC IMMUNITY

Immunity Test	IEC 60601 Test level	Compliance	Electromagnetic Environment-Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be greater than 30%.
Electrical fast transient/burst IEC61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or residential environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Not applicable	Mains power quality should be that of a typical commercial or residential environment.
Voltage dips, short interruptions, and voltage variations on power supply lines IEC 61000-4-11	< 5% <i>UT</i> (>95% dip in <i>UT</i> ) for 0.5 cycle  40% <i>UT</i> (60% dip in <i>UT</i> ) for 5 cycles  70% <i>UT</i> (30% dip in <i>UT</i> ) for 25 cycles  <5% <i>UT</i> (>95% dip in <i>UT</i> ) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or residential environment. If the user of the system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptable power supply or its battery.
Power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or residential environment.

**Note:** *UT* is the A.C. mains voltage prior to application of the test level.

Immunity test	IEC / EN 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC / EN 61000-4-6</p> <p>Radiated RF IEC / EN 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any parts of the system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>d = 1,2\sqrt{P}</math></p> <p><math>d = 1,2\sqrt{P}</math> 80 MHz to 800 MHz</p> <p><math>d = 2,3\sqrt{P}</math> 800 MHz to 2,5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (a) should be less than the compliance level in each frequency range (b)</p>

			<p>Interference may occur near equipment marked with the following symbol:</p> 
<p><b>NOTE 1</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			

<sup>(a)</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the system.

<sup>(b)</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**NOTICE:** There are no deviations from the collateral standard and allowances uses

**NOTICE:** All necessary instruction for maintaining compliance regarding EMC can be found in the general maintenance section in this instruction. No further steps required.

Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. To ensure compliance with the EMC requirements as specified in IEC 60601-1-2, it is essential to use only the following accessories.

Item	Manufacturer	Item #
IP30 Insert Earphone	RadioEar	8101884
DD45 Headphones	RadioEar	8503851
DD450 Headphones	RadioEar	8520321
DD65 v2 Headphones	RadioEar	8510115
B81 Bone Vibrator	RadioEar	8501119

Conformance to the EMC requirements as specified in IEC 60601-1-2 is ensured if the cable types and cable lengths are as specified below:

Description	Length	Screened/Unscreened
IP30 Insert Earphone	2.5 m	Unscreened
DD45 Headphones	2.5 m	Unscreened
DD450 Headphones	2.5 m	Unscreened
DD65 v2 Headphones	2.5 m	Unscreened
B81 Bone Vibrator	2 m	Unscreened
Power Supply Cable	2.5 m	Unscreened
USB Cable	2 m	Screened

**NOTICE:** The use of the accessories, transducers, and cables with medical equipment/system other than this equipment may result in increased emissions or decreased immunity of the medical equipment/system.

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The screenshot displays the GSI AudioStar Pro software interface for a speech test. At the top, it shows 'Patient ...', 'Speech Table', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 shows a score of 50 dB HL, and Channel 2 shows a score of 0 dB HL. A central 'Playground' area is visible. Below the channels are sections for 'Speech Test Results' and a list of word lists. The word list table has columns for Ear, Test Type, Int Ext Mic, Word Lists, Aid, %, dB HL, and dB EM. The word list table shows various words like Daybreak, Northwest, Mushroom, and Doormat. The interface also includes a 'Word Nav' button and a 'dB Step' of 5. The bottom right corner shows the time 2:05 PM and date 02-25-2019.

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
## BKB-SIN

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' with 'R', 'B', and 'L' indicators. There are also fields for 'INT A' and 'INT B' with 'Right' and 'Left' labels. Below these are 'PTA AC: BC' and 'SII' fields for both sides. The 'Test Results Group 1' section contains three tables for 'Right', 'Left', and 'Reliability' ears, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. The 'Word List' section shows a list of 10 word pairs (A1-A10) with their corresponding scores (S/N). The bottom status bar includes 'Word Lists', 'Word Nav', 'Age', 'Aided', 'dB Step' (5), 'Group' (1), 'Research', and a timestamp '10:03 AM 02-26-2019'.

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## QUICKSIN

The screenshot shows the QuickSIN software interface. At the top, there are two channels: Channel 1 (70 dB HL) and Channel 2 (0 dB HL). Each channel has a section for SNR Loss Averages with tables for Basic, HFE, and HFE-LP tests for Right (R), Back (B), and Left (L) ears. Below this are Test Results Group 1 tables for both ears, with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. A word list section is visible, showing a track of 3 with the following text: "A WHITE SILK JACKET goes with ANY SHOES", "The CHILD CRAWLED INTO the DENSE GRASS", "FOOTPRINTS SHOWED the PATH he TOOK up the BEACH", "A VENT NEAR the EDGE brought in FRESH AIR", "It is a BAND of STEEL THREE INCHES WIDE", and "The WEIGHT of the PACKAGE was SEEN on the HIGH SCALE". A score table on the right shows S/N ratios from 25 to 0, all with a score of "-", and a Sum row. The bottom bar contains various settings: Word Lists, Word Nav, Aided, dB Step (5), Group (1), Research, and a timestamp of 10:53 AM on 03-01-2019.

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### Aktualizace

The screenshot shows the 'Configuration' screen with the following data:

Instrument			
App Version	2.0.0 (Build 689)	Login	No Login Required
DSPAud Version	2.21 (04/01/2019 12:47)	Language	English (United States)
OS Version	2018, 1, 2018	Symbol Set	USA
Serial Number	PP10266	Bone	Mastoid
License number:	48536688	Assistant Headset	On
Free Disk Space	78.84 MB	Display Brightness	85
Free SD Card Space	9.89 GB	Mode	AWL
Last Calibration	01/17/2019	Next Calibration	12/10/2019

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Update | Bone: Mastoid | Assistant Headset:  | Export Logs | 3:50 PM 05/28/2019

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## BKB-SIN

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into two channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Below the channels are sections for SNR Loss Averages (Group 1 and Group 2) with R, B, and L indicators. There are also fields for PTA AC: BC, SII, and Reliability for both Right and Left ears. The interface includes three 'Test Results Group 1' tables with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. A central 'BKB-SIN : List Pair 1' table contains 10 word pairs (A1-A10). On the sides, there are 'Score' tables with columns for S/N and a dash. The bottom bar contains controls for Word Lists, Word Nav, Age, Aided, dB Step (set to 5), Group (set to 1), Research, and a timestamp of 10:03 AM on 02-26-2019.

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## QUICKSIN

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (red header) and Channel 2 (blue header). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' section with tables for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below the channels are 'Test Results Group 1' tables for each ear, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains several sentences for testing. To the right of the list is a 'Score' table with columns for 'S/N' and 'Score'. At the bottom, there is a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

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The screenshot displays the AudioStar Pro software interface for a speech test. At the top, it shows 'Patient ...', 'Speech Table', and 'Examiner 1'. The main area is divided into 'Channel 1' (red) and 'Channel 2' (blue). Channel 1 shows a score of 50 dB HL, and Channel 2 shows 0 dB HL. A 'Score - SRT' section indicates 0%. Below this is a 'Playground' section. The interface includes fields for 'PTA AC: BC', 'SII', and 'Reliability'. Two 'Speech Test Results' tables are present, one for each channel, with columns for Ear, Test Type, Int Ext Mic, Word Lists, Aid, %, dB HL, and dB EM. At the bottom, a 'Word Lists' table is shown with words like Daybreak, Northwest, Mushroom, and Doormat. The interface also displays 'PTA AC: BC', 'SII', and 'Reliability' fields. The bottom status bar shows 'Test Type', 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (5), and the time '2:05 PM 02-25-2019'.

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
## BKB-SIN

The screenshot displays the BKB-SIN software interface. At the top, it shows 'Patient ...', 'BKB-SIN', and 'Examiner 1'. The main area is divided into 'Channel 1' (70 dB HL) and 'Channel 2' (0 dB HL). Below the channels are 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages' for Right (R), Back (B), and Left (L) ears. There are also 'INT A' and 'INT B' sections for 'Right' and 'Left' ears, each with a 'Phone' option and a scale from -20 to 3. Below these are 'PTA AC: BC' and 'SII' fields for both ears, along with a 'Reliability' field. The 'Test Results Group 1' section contains three empty tables with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. At the bottom, there is a 'BKB-SIN : List Pair 1' section with a table of word pairs and their scores. The interface also includes a 'Word Lists' section and a status bar at the bottom with various settings and a timestamp of 10:03 AM on 02-26-2019.

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## QUICKSIN

The screenshot displays the QuickSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. Below this, there are two main channels: Channel 1 (Right) and Channel 2 (Left). Channel 1 is set to 70 dB HL, and Channel 2 is set to 0 dB HL. Each channel has a frequency scale from -20 to 3 dB HL. The interface also includes sections for 'Group 1 SNR Loss Averages' and 'Group 2 SNR Loss Averages', each with a table for 'Basic', 'HFE', and 'HFE-LP' across 'R', 'B', and 'L' ears. Below these are 'Test Results Group 1' tables for each ear, with columns for 'Ear', 'Word List', 'dB HL', 'SNR 50', and 'SNR Loss'. A 'QuickSin : List 1 (Track 3)' section contains a list of sentences for testing. To the right of this list is a 'Score' table with columns for 'S/N' and 'Score'. At the bottom, there is a control bar with buttons for 'Word Lists', 'Word Nav', 'Aided', 'dB Step' (set to 5), 'Group' (set to 1), 'Research', and a clock showing '10:53 AM 03-01-2019'.

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The screenshot shows the 'Configuration' screen for a patient. The top bar includes 'Patient ...' and 'Configuration'. Below this is a table titled 'Instrument' with the following data:

Instrument			
App Version	2.0.0 (Build 689)	Login	No Login Required
DSPAud Version	2.21 (04/01/2019 12:47)	Language	English (United States)
OS Version	2018, 1, 2018	Symbol Set	USA
Serial Number	PP10266	Bone	Mastoid
License number:	48536688	Assistant Headset	On
Free Disk Space	78.84 MB	Display Brightness	85
Free SD Card Space	9.89 GB	Mode	AWL
Last Calibration	01/17/2019	Next Calibration	12/10/2019

Below the instrument table is another table for 'Transducer' settings:

Transducer		Right	Left
Headset Phone	DD45	45r	45l
Insert Phone	IP30	ip30r	ip30l
High Frequency Phone	DD450	450r	450l
Bone Conductor	B81	81	

At the bottom of the screen, there is a footer with contact information for 'Audiology of Illinois' and a status bar with 'Update', 'Bone Mastoid', 'Assistant Headset', 'Export Logs', a key icon, a gear icon, and the time '3:50 PM 05/28/2019'.

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The screenshot shows the GSI AudioStar Pro software interface. At the top, there are fields for 'Patient ...', 'Speech Table', and 'Examiner 1'. Below this, there are two main channels: 'Channel 1' (red) and 'Channel 2' (blue). Channel 1 shows a score of 50 dB HL, and Channel 2 shows a score of 0 dB HL. A central 'Playground' area is visible. Below the channels are 'Speech Test Results' tables for both channels. At the bottom, there is a 'Basic Auditory Tests - Adult : Spondee A' list with words like Daybreak, Northwest, Mushroom, and Doormat. The interface also shows a 'Test Type' dropdown set to 'Word Lists', 'Word Nav', 'Aided', and 'dB Step' set to 5. The bottom right corner shows the time 2:05 PM and date 02-25-2019.

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
**Channel 2** 0 dB HL  
INT B Left Phone

PTA AC: BC: SII: Right Reliability: PTA AC: BC: SII: Left

Test Results Group 1					Test Results Group 1					Test Results Group 1				
Ear	Word List	dB HL	SNR 50	SNR Loss	Ear	Word List	dB HL	SNR 50	SNR Loss	Ear	Word List	dB HL	SNR 50	SNR Loss

Score	BKB-SIN : List Pair 1		Page 1/2	Score	
S/N 21 -	A1	THEY are LOOKING AT the CLOCK	A6	HE PLAYED with his TRAIN	S/N 6 -
S/N 18 -	A2	The CAR ENGINE is RUNNING	A7	The BAG FELL to the GROUND	S/N 3 -
S/N 15 -	A3	CHILDREN LIKE STRAWBERRIES	A8	The BOY DID a HANDSTAND	S/N 0 -
S/N 12 -	A4	THEY are BUYING some BREAD	A9	The WATER BOILED QUICKLY	S/N -3 -
S/N 9 -	A5	The GREEN TOMATOES are SMALL	A10	The MAN is PAINTING a SIGN	S/N -6 -
					Sum -

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## QUICKSIN

The screenshot displays the QUICKSIN software interface. At the top, it shows 'Patient ...' and 'QuickSIN'. The main area is split into two channels: Channel 1 (70 dB HL) and Channel 2 (0 dB HL). Each channel has a frequency scale from -20 to 3 dB HL. Below the channels are sections for SNR Loss Averages (Group 1 and Group 2) with tables for Basic, HFE, and HFE-LP tests for Right (R), Both (B), and Left (L) ears. There are also PTA AC: BC and SII: fields for both Right and Left ears. The Test Results Group 1 section contains three tables with columns for Ear, Word List, dB HL, SNR 50, and SNR Loss. Below these is a 'QuickSin : List 1 (Track 3)' section with a list of sentences. To the right is a 'Score' table with columns for S/N and a dash. The bottom bar contains various controls: Word Lists, Word Nav, Aided, dB Step (5), Group (1), Research, and a timestamp of 10:53 AM on 03-01-2019.

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Basic: Auditory Tests - Adult : Spondee A Page 1/2

Word Lists	Daybreak	Northwest	Mushroom	Doormat
Playground	Iceberg	Padlock	Sunset	Duck pond
Eardrum	Inkwell	Baseball	Whitewash	Oatmeal
Cowboy	Hot dog	Mousetrap	Airplane	Headlight
Greyhound	Stairway	Woodwork	Drawbridge	Armchair
Hothouse				

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
Score BKB-SIN : List Pair 1 Page 1/2

Score	Word List	Score
S/N 21 -	A1 THEY are LOOKING AT the CLOCK	A6 HE PLAYED with his TRAIN
S/N 18 -	A2 The CAR ENGINE is RUNNING	A7 The BAG FELL to the GROUND
S/N 15 -	A3 CHILDREN LIKE STRAWBERRIES	A8 The BOY DID a HANDSTAND
S/N 12 -	A4 THEY are BUYING some BREAD	A9 The WATER BOILED QUICKLY
S/N 9 -	A5 The GREEN TOMATOES are SMALL	A10 The MAN is PAINTING a SIGN
S/N 6 -		
S/N 3 -		
S/N 0 -		
S/N -3 -		
S/N -6 -		
Sum -		

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QuickSin : List 1 (Track 3)

A WHITE SILK JACKET goes with ANY SHOES
The CHILD CRAWLED INTO the DENSE GRASS
FOOTPRINTS SHOWED the PATH he TOOK up the BEACH
A VENT NEAR the EDGE brought in FRESH AIR
It is a BAND of STEEL THREE INCHES WIDE
The WEIGHT of the PACKAGE was SEEN on the HIGH SCALE

Score	
S/N 25	-
S/N 20	-
S/N 15	-
S/N 10	-
S/N 5	-
S/N 0	-
Sum	-

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**Copyright © 2021 Grason-Stadler.** Tous droits réservés. Aucune partie de cette publication ne peut être reproduite ou transmise sous quelque forme ou par quelque moyen que ce soit sans l'autorisation écrite préalable de Grason-Stadler. Les informations contenues dans cette publication sont la propriété de Grason-Stadler.
**Copyright © 2021 Grason-Stadler.** Minden jog fenntartva. A kiadvány semmilyen formában és eszközzel nem sokszorosítható vagy továbbítható a Grason-Stadler előzetes írásbeli engedélye nélkül. A jelen kiadványban szereplő információk a Grason-Stadler tulajdonát képezik.
**Copyright © 2021 Grason-Stadler.** Tutti i diritti riservati. Nessuna parte di questa pubblicazione può essere riprodotta o trasmessa in alcuna forma o con alcun mezzo senza la previa autorizzazione scritta di Grason-Stadler. Le informazioni contenute nella presente pubblicazione sono di proprietà di Grason-Stadler.