

Instructions for use

IIC-CIC-ITC
ITE HS-ITE FS

Oticon Opn
Oticon Siya



Made for

iPhone | iPad | iPod

oticon
life-changing technology



WARNING: People younger than 18 should go to a doctor before using this.

People younger than 18 years old need specialized care, and using this without a medical evaluation may worsen impairment or disability. A hearing aid user who is younger than 18 should have a recent medical evaluation from a doctor, preferably an ear-nose-throat doctor (an ENT). Before using this, a doctor should determine that the use of a hearing aid is appropriate.

WARNING to Hearing Aid Dispensers:

You should advise a prospective hearing aid user to consult promptly with a doctor, preferably an ear specialist such as an ENT, before dispensing a hearing aid if you determine through inquiry, actual observation, or review of any other

available information concerning the prospective user, that the prospective user has any of the following conditions:

- Visible deformity of the ear, either congenital or traumatic
- Fluid, pus, or blood coming out of the ear within the previous 6 months
- Pain or discomfort in the ear
- History of excessive ear wax or suspicion that something is in the ear canal
- Dizziness, either recent or long-standing
- Sudden, quickly worsening, or fluctuating hearing loss within the previous 6 months
- Hearing loss or ringing (tinnitus) only in one ear or a noticeable difference in hearing between ears
- Audiometric air-bone gap equal to or greater than 15 dB at 500 Hz, 1000 Hz, and 2000 Hz

**WARNING to Hearing Aid Dispenser,
Outputs over 132 dB SPL:**

You should exercise special care in selecting and fitting a hearing aid with a maximum output that exceeds 132 dB SPL because it may impair the remaining hearing of the hearing aid user.

Caution: This is not hearing protection.

You should remove this device if you experience overly loud sounds, whether short or long-lasting. If you're in a loud place, you should use the right kind of hearing protection instead of wearing this device. In general, if you would use ear plugs in a loud place, you should remove this device and use ear plugs.

Caution: The sound output should not be uncomfortable or painful.

You should turn down the volume or remove the device if the sound output is uncomfortably loud or painful. If you consistently need to turn the volume down, you may need to further adjust your device.

Caution: You might need medical help if a piece gets stuck in your ear.

If any part of your hearing aid, like the eartip, gets stuck in your ear, and you can't easily remove it with your fingers, get medical help as soon as you can. You should not try to use tweezers or cotton swabs because they can push the part farther into your ear, injuring your eardrum or ear canal, possibly seriously.

Note: What you might expect when you start using a hearing aid

A hearing aid can benefit many people with hearing loss. However, you should know it will not restore normal hearing, and you may still have some difficulty hearing over noise. Further, a hearing aid will not prevent or improve a medical condition that causes hearing loss.

People who start using hearing aids sometimes need a few weeks to get used to them. Similarly, many people find that training or counseling can help them get more out of their devices.



If you have hearing loss in both ears, you might get more out of using hearing aids in both, especially in situations that make you tired from listening –for example, noisy environments.

Note: Tell FDA about Injuries, malfunctions, or other adverse events.

To report a problem involving your hearing aid, you should submit Information to FDA as soon as possible after the problem. FDA calls them “adverse events,” and they might include: skin irritation in your ear, injury from the device (like cuts or scratches, or burns from an overheated battery), pieces of the device getting stuck in your ear, suddenly worsening hearing loss from using the device, etc.

Instructions for reporting are available at <https://www.fda.gov/Safety/MedWatch>, or call 1-800-FDA-1088. You can also download a form to mail to FDA.

Note: Hearing loss in people younger than 18

- People younger than 18 should see a doctor first, preferably an ear-nose-throat doctor (an ENT), because they may have different needs than adults.
- The doctor will identify and treat medical conditions as appropriate.
- The doctor may refer the person to an audiologist for a separate test, a hearing aid evaluation.
- The hearing aid evaluation will help the audiologist select and fit the appropriate hearing aid.

A person who is younger than 18 years old with hearing loss should have a medical evaluation by a doctor, preferably an ENT, before buying a hearing aid. The purpose of a medical evaluation is to identify and treat medical conditions that may affect hearing but that a hearing aid won't treat on its own.

Following the medical evaluation and if appropriate, the doctor will provide a written statement that the hearing loss has been medically evaluated and the person is a candidate for a hearing aid. The doctor may refer the person to an audiologist for a hearing aid evaluation, which is different from the medical evaluation and is intended to identify the appropriate hearing aid.

The audiologist will conduct a hearing aid evaluation to assess the person's ability to hear with and without a hearing aid. This will enable the audiologist to select and fit a hearing aid for the person's individual needs. An audiologist can also provide evaluation and rehabilitation since, for people younger than 18, hearing loss may cause problems in language development and educational and social growth. An audiologist is qualified by training and experience to assist in the evaluation and rehabilitation of hearing loss in people younger than 18.

Model overview

This booklet is valid for the Oticon Opn™ and Oticon Siya families in the following hearing aid models, battery sizes and styles:

Oticon Opn FW6:

- Oticon Opn 1 GTIN: (01) 05707131340795
- Oticon Opn 2 GTIN: (01) 05707131340818
- Oticon Opn 3 GTIN: (01) 05707131340801

Oticon Siya FW1:

- Oticon Siya 1 GTIN: (01) 05707131340887
- Oticon Siya 2 GTIN: (01) 05707131340894

Battery size:

10 312 13

Bluetooth®: Yes No

Styles:

- IIC Invisible-In-the-Canal
- CIC Completely-In-the-Canal
- ITC In-the-Canal
- ITE HS In-the-Ear Half Shell
- ITE FS In-the-Ear Full Shell

Introduction to this booklet

This booklet gives you guidance on how to use and maintain your new hearing aid. Please read the booklet carefully including the **Warning section**. This will help you to achieve the full benefit of your new hearing aid.

Your hearing care professional has adjusted the hearing aid to meet your needs. If you have additional questions, please contact your hearing care professional.

| **About** | Startup | Handling | Options | Tinnitus | Warnings | More info |

For your convenience this booklet contains a navigation bar to help you navigate easily through the different sections.

Intended use

The hearing aid is intended to amplify and transmit sound to the ear and thereby compensate for impaired hearing within mild to severe hearing loss. This hearing aid is intended for use by adults and children older than 36 months.

IMPORTANT NOTICE

The hearing aid amplification is uniquely adjusted and optimized to your personal hearing capabilities during the hearing aid fitting performed by your hearing care professional.

Table of contents

About

Identify your hearing aid	12
Identify left and right hearing aid	16
MultiTool for handling batteries and cleaning	17

Start up

Turn the hearing aid on and off	18
Replace the battery	19

Handling

Insert the hearing aid	22
Caring for your hearing aid	24
Filter replacement	26
Flight mode (hearing aids with Bluetooth and push-button)	30
Optional features and accessories	31

Options

Change programs	32
Change volume with push button	34
Change volume with volume wheel	35
Mute the hearing aid	36


	Use hearing aids with iPhone, iPad, and iPod touch	37
	Wireless accessories	38
	Other options	40
Tinnitus		
	Tinnitus SoundSupport™ (optional)	42
	Guidelines for tinnitus sound generator users	43
	Sound options and volume adjustment	44
	Limitation on use time	49
	Warnings related to Tinnitus	52
Warnings		
	General warnings	53
More info		
	Troubleshooting guide	58
	Cell phone	60
	Water & dust resistant (IP68)	62
	Conditions of use	63
	Warranty	64
	Technical information	66
	Your individual hearing aid settings	72
	Summary of relevant studies	74
	Technical Data	76

Identify your hearing aid

For your in-the-ear hearing aid, there are several different styles available with different battery sizes. Please see the model overview for your style and battery size. This will make it easier for you to navigate through this booklet.

Size 10 battery (CIC shown)


What it is and does



Components may be positioned differently on your hearing aid.

Size 312 battery (ITC shown)

What it is and does




Size 13 battery (ITE HS shown)

What it is and does

Push button (optional)
Mute, change volume and program


Volume wheel (optional)
Change volume



Microphone inlet
Sound in

Vent

Battery drawer
Contains the battery and functions as on/off switch




Vent

Sound outlet
Sound out

Components may be positioned differently on your hearing aid.

Identify left and right hearing aid


It is important to distinguish between the left and the right hearing aid as they might be shaped and programmed differently.



* Shell color and dot is only available for IIC

MultiTool for handling batteries and cleaning

The MultiTool contains a magnet that makes it easier to replace the battery in the hearing aid. It also contains a brush and a wire loop for cleaning and removing earwax. If the vent is very small, a special tool may be required. If you need a new tool please contact your hearing care professional.



IMPORTANT NOTICE


The MultiTool has a built-in magnet. Keep the MultiTool 30 cm away from credit cards and other magnetically sensitive devices.

Turn the hearing aid on and off

The battery drawer is also used to switch the hearing aid on and off. To save battery life, make sure your hearing aid is switched off when you are not wearing it. If you wish to return to the standard settings of the hearing aid, simply open and then close the battery drawer (quick reset).


Turn ON

Close the battery drawer with the battery in place.



Turn OFF

Open the battery drawer.



Note: Open the battery door fully to make sure the hearing aid does not drain the battery fully, and to allow air to circulate whenever you are not using your hearing instrument, especially at night or for longer periods of time!

Replace the battery

When it is time to replace the battery, you will hear three beeps repeated at moderate intervals until the battery runs out.



Three beeps*

= The battery is running low



Four beeps

= The battery has run out

Battery tip


To make sure the hearing aid is always working, bring spare batteries with you, or replace the battery before you leave home.

* If your hearing aid has Bluetooth® this will be turned off and it will not be possible to use wireless accessories.

Note: Batteries need to be replaced more often if you are streaming audio or music to your hearing aids.

How to replace the battery

1. Remove



Fully open the battery drawer. Remove the battery.


2. Uncover



Remove the sticky label from the + side of the new battery.


Tip:
Wait 2 minutes so that the battery can draw air, to ensure optimal functioning.

3. Insert




Insert the new battery into the battery drawer. Make sure the + side is facing up.

4. Close



Close the battery drawer. The hearing aid will play a jingle. Hold the hearing aid close to your ear to hear the jingle.

Tip



The MultiTool can be used for battery change. Use the magnetic end to remove and insert batteries.

The MultiTool is provided by your hearing care professional.


Insert the hearing aid

Step 1



Place the tip of the hearing aid in your ear canal.

Step 2



Gently pull your ear outwards and push the hearing aid into the ear canal, twisting slightly if necessary. Follow the natural contour of the ear canal. Push the hearing aid to make sure it fits comfortably in the ear.

Remove your hearing aid

Hold the hearing aid by the pull-out string (if available). Gently pull the hearing aid from the ear canal.

If your hearing aid doesn't have a pull-out string, you can remove it by pulling on the edge of the hearing aid.

IMPORTANT NOTICE


DO NOT use the battery door as a handle to insert or remove your hearing aid. It is not designed for this purpose.

Caring for your hearing aid



When handling your hearing aid, hold it over a soft surface to avoid damage if you drop it.

Cleaning the hearing aid

Carefully brush away debris from the microphone inlets with a clean brush. Gently brush the surface. Make sure that the filters do not fall off.



Clean the vent by pressing the brush through the hole while twisting it slightly.



IMPORTANT NOTICE

Use a soft, dry cloth to clean the hearing aid. It must never be washed or immersed in water or other liquids.

Filter replacement

The filters keep wax and debris from damaging the hearing aid. If the filters become clogged, please replace the filters or contact your hearing care professional.

- Prowax miniFit filter protects the sound outlet.
- O-cap and T-cap protect the microphone inlet.

Please refer to the following pages for instructions on how to replace the appropriate filters.


IMPORTANT NOTICE

Always use the same type of wax filter as was originally supplied with the hearing aid.

If you are in any doubt about the use or replacement of wax filters, contact your hearing care professional.


□ Replace ProWax miniFit filter

1. Tool




Remove the tool from the shell. The tool has two pins, one empty for removal and one with the new ProWax miniFit filter.

2. Remove



Push the empty pin into the ProWax miniFit filter in the hearing aid and pull it out.


3. Insert



Insert the new ProWax miniFit filter using the other pin, remove the tool and throw it out.


□ **Replace T-Cap filter** (hearing aids with 10 batteries)

1. Tool




Remove the tool from the packaging. The tool has two ends, one for removal and one with the new T-Cap filter.

2. Remove



Push the tool fork under the top edge of the used T-Cap filter and lift it out.


3. Insert



Insert the new T-Cap filter and remove the tool by twisting it slightly. Throw the tool out after use.


□ **Replace O-Cap filter** (hearing aids with 312 and 13 batteries)

1. Tool




Remove the tool from the packaging. The tool has two ends, one for removal and one with the new O-Cap filter.

2. Remove



Push the pointed end of the tool into the existing O-Cap filter and pull it out.

3. Insert



Insert the new O-Cap filter using the other end of the tool, remove the tool and throw it out.

□ **Flight mode**


(hearing aids with Bluetooth and push-button)

When boarding an airplane or entering an area in which it is prohibited to radiate radio signals, e.g. during flight, flight mode must be activated. The hearing aid will still be working. It is only necessary to activate flight mode on one hearing aid, to turn off Bluetooth on both hearing aids. If your hearing aid does not have a push-button you need to turn off your hearing aid.

To activate and deactivate

Press the push button for at least 7 seconds. A jingle confirms your action.

Opening and closing the battery drawer will also deactivate flight mode.



Optional features and accessories


The features and accessories described on the following pages are optional. Please contact your hearing care professional to find out how your hearing aid is programmed.

If you experience difficult listening situations, a special program may be helpful. These are programmed by your hearing care professional.

Write down any hearing situations in which you may need help.

□ **Change programs**

Your hearing aid can have up to 4 different programs.
These are programmed by your hearing care professional.







• Press the button to change program. Use a short press if the push button is used for program change only and a long press if it is also used for volume control.

Note that if you have two synchronized hearing aids, (both hearing aids respond when either push button is operated) the RIGHT hearing aid switches forward from e.g., program 1 to 2 and the LEFT hearing aid switches backwards from e.g., program 4 to 3.

If your hearing aids work independently, you must press the buttons on each hearing aid.


To be filled out by the hearing care professional

Program	Sound you will hear when activated	When to use
1	 "1 beep"	
2	 "2 beeps"	
3	 "3 beeps"	
4	 "4 beeps"	

Program change:	
<input type="checkbox"/> Independent	<input type="checkbox"/> Synchronized
<input type="checkbox"/> LEFT	<input type="checkbox"/> RIGHT
<input type="checkbox"/> Short press	<input type="checkbox"/> Long press

Change volume with push button

The push button allows you to adjust the volume.
You may hear a click when you turn the volume up or down.



A short press on the
RIGHT hearing aid
increases the volume

A short press on the
LEFT hearing aid
decreases the volume



To be filled out by the hearing care professional


Volume change

LEFT

RIGHT


□ **Change volume with volume wheel**

The volume wheel allows you to adjust the volume.
You may hear a click when you turn the volume up or down.




Turn forward to
turn up volume

Turn backwards
to turn down the
volume



Mute the hearing aid

Use the mute function if you need to silence the hearing aid. Only available for hearing aids with push button.



Apply a very long press to the push button to mute the hearing aid. To reactivate the hearing aid, push the button briefly.

NOTE: the mute function only mutes the microphone(s) on the hearing aids.

IMPORTANT NOTICE

Do not use the mute function as an off switch, as the hearing aid still draws current from the battery in this mode.

Use hearing aids with iPhone, iPad, and iPod touch

Your hearing aids are Made for iPhone® and allow for direct communication and control with iPhone, iPad® or iPod touch®. For assistance in using these products with your hearing aids,

please contact your hearing care professional or visit:
www.oticon.com/connectivity



To see how to pair your hearing aid with your iPhone visit
www.oticon.com/support

For information on compatibility, visit www.oticon.com/support/compatibility

Use of the Made for Apple badge means that an accessory has been designed to connect specifically to the Apple products identified in the badge, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards.

Note that use of this Accessory with iPod, iPhone, or iPad may affect wireless performance.

Wireless accessories

As an enhancement to your hearing aids, a range of wireless accessories are available. These can enable you to hear and communicate better in many everyday situations.

ConnectClip

When paired with your mobile phone, you can use the hearing aids as a hands-free headset, or as a remote microphone.

TV Adapter 3.0

Wireless transmitter of sound from TV and electronic audio devices. TV Adapter streams sound directly to your hearing aids.

Remote Control 3.0

Offers the ability to change programs, adjust volume, or mute your hearing aids.

Oticon ON App

Offers an intuitive and discreet way to control your hearing aids. For iPhone, iPad, iPod touch, and Android™ devices.

□ **Phone Adapter 2.0**

Connects wirelessly to the hearing aids via ConnectClip allowing for use of landline phones.

For more information visit: www.oticon.com/connectivity

Other options



Telecoil - only optional for some styles*

Telecoil helps you hear better when using a telephone with a built-in loop or when you are in buildings with teleloop systems such as theaters, churches, or lecture rooms. This symbol or a similar sign is shown wherever a teleloop has been installed.



Autophone - only optional for some styles

The autophone can automatically activate a phone program in the hearing aid, if you place a dedicated magnet on your phone.

For more information, please contact your hearing care professional.

Tinnitus SoundSupport™ (optional)

Intended use of Tinnitus SoundSupport

Tinnitus SoundSupport is a tool intended to generate sounds to provide temporary relief for patients suffering from tinnitus as part of a tinnitus management program. The target population is the adult population (over 18 years old).

Tinnitus SoundSupport is targeted to licensed hearing care professionals (audiologists, hearing aid specialists, or otolaryngologists) who are familiar with the evaluation and treatment of tinnitus and hearing loss. Fitting of Tinnitus SoundSupport must be done by a hearing care professional participating in a tinnitus management program.

Styles available with Tinnitus SoundSupport

CIC, ITC, ITE HS & ITE FS for hearing aids with push button.

Please see the model overview to find out which style your hearing aid is.

Guidelines for tinnitus sound generator users

These instructions contain information about Tinnitus SoundSupport, which may have been enabled in your hearing aids by your hearing care professional.

Tinnitus SoundSupport is a tinnitus management tool intended to generate sound of sufficient intensity and bandwidth to help manage tinnitus.

Your hearing care professional will also be able to offer the appropriate follow-up care. It is important to follow his/her advice and directions regarding such care.

Prescription use only

Good health practice requires that a person reporting tinnitus has a medical evaluation by a licensed ear physician before using a sound generator. The purpose of such an evaluation is to ensure that any medically treatable condition that may cause tinnitus is identified and treated prior to using a sound generator.

Sound options and volume adjustment

Tinnitus SoundSupport is programmed by your hearing care professional to match your hearing loss and preferences for tinnitus relief. It offers a number of different sound options. Together with your hearing care professional, you can select your preferred sound(s).

Tinnitus SoundSupport programs

Together with your hearing care professional you decide for which programs you may want to have Tinnitus SoundSupport activated. The sound generator can be activated in up to four different programs.

Mute

If you are in a program for which Tinnitus SoundSupport is activated, the mute functionality will mute only the environmental sounds, and not the sound from Tinnitus SoundSupport.

See section: "Mute the hearing aid".

Volume adjustments with Tinnitus SoundSupport

When Tinnitus SoundSupport is activated (in a program), you will only be able to adjust the volume of Tinnitus SoundSupport using the push button or volume wheel. This means that you cannot adjust the environmental sounds in this program.

Your hearing care professional will set the volume control for Tinnitus SoundSupport in one of two ways:

- A) Change volume in each ear separately, or
- B) Change volume in both ears simultaneously.

Change volume with push button

See section “Change volume with push button” for illustration

A) How to change Tinnitus SoundSupport volume in each ear separately

To **increase** volume (on one hearing aid only), use a short press on the push button repeatedly until desired level is reached. The sound will always be louder with the first press(es) until two beeps are heard. Hereafter the volume will decrease.

To **decrease** volume (on only one hearing aid), continue to press the push button repeatedly until desired level is reached.

B) How to change Tinnitus SoundSupport volume in both ears simultaneously

You can use one hearing aid to increase the sound and the other hearing aid to decrease the sound:

To **increase** volume, use a short press on the push button repeatedly on the RIGHT hearing aid.

To **decrease** volume, use a short press on the push button repeatedly on the LEFT hearing aid.

To be filled out by your hearing care professional.

Change volume with volume wheel

See section “Change volume with volume wheel” for illustration

A) How to change Tinnitus SoundSupport volume in each ear separately

To **increase** volume (on one hearing aid only), turn the volume wheel forward.

To **decrease** volume (on one hearing aid only), turn the volume wheel backwards.

B) How to change Tinnitus SoundSupport volume in both ears simultaneously

You can use one hearing aid to increase/decrease the sound in both hearing aids. When changing the volume in one hearing aid, the volume on the other hearing aid will follow.

To **increase** volume, turn the volume wheel forward.

To **decrease** volume, turn the volume wheel backwards.

To be filled out by your hearing care professional.

Hearing aid with Bluetooth

If your hearing aid has Bluetooth connection, you will be able to adjust the environmental sounds in a program with Tinnitus SoundSupport with Remote Control 3.0 or ConnectClip.

Oticon ON App for iPhone, iPad, iPod touch, and Android™ devices offers the possibility to adjust both the environmental sound and the Tinnitus SoundSupport in each program directly from the app.

Limitation on use time

Daily use

The volume levels of Tinnitus SoundSupport can be set to a level which could lead to permanent hearing damage when used for a prolonged period of time. Your hearing care professional will advise you of the maximum amount of time per day you should use Tinnitus SoundSupport. It should never be used at uncomfortable levels.

See table “Tinnitus SoundSupport: Limitation on use” in the section “Your individual hearing aid settings” at the end of this booklet to learn how many hours per day you can safely use the relief sound in your hearing aids.

Important information for hearing care professionals about Tinnitus SoundSupport

Device description

Tinnitus SoundSupport is a module function that can be enabled in the hearing aids by the hearing care professional.

Maximum wearing time

The wearing time of Tinnitus SoundSupport will decrease as you increase the level above 80 dB(A) SPL. The fitting software will automatically display a warning when the hearing aid exceeds 80 dB(A) SPL. See “Max wearing time indicator” next to the tinnitus fitting graph in the fitting software.

The volume control is deactivated

By default the volume control for the sound generator is deactivated in the hearing aid. Risk of noise exposure increases when the volume control is activated.

If the volume control is activated

A warning may be displayed if you activate the tinnitus volume control in the “Buttons and Indicators” screen. This occurs if the relief sound can be listened to at levels that may cause hearing damage. The “Max wearing time” table in the fitting software displays the number of hours the patient can safely use Tinnitus SoundSupport.

- Note the max wearing time for each program for which Tinnitus SoundSupport is activated.
- Write those values in the table: “Tinnitus SoundSupport: Limitation on use”, in the back of this booklet.
- Instruct your patient accordingly.

Warnings related to Tinnitus

If your hearing care professional has activated the sound generator Tinnitus SoundSupport, please pay attention to the following warnings.

There are some potential concerns associated with the use of any sound generated by a tinnitus management device. Among them are the potential worsening of tinnitus, and/or a possible change in hearing thresholds.

Should you experience or notice a change in hearing or tinnitus, or any dizziness, nausea, headaches, heart palpitations, or possible skin irritation at the point of contact with the device, you should immediately discontinue use of the device and consult a medical, audiology, or other hearing care professional.

As with any device, misuse of the sound generator feature may cause potentially harmful effects. Care should be taken to prevent unauthorized use and to keep the device out of reach of children and pets.

Maximum wearing time

Always follow the maximum wearing time per day of the Tinnitus SoundSupport advised by your hearing care professional. Prolonged use may lead to worsening of your tinnitus or of your hearing loss.

General warnings

For your personal safety and to ensure correct usage, you should familiarize yourself fully with the following general warnings before using your hearing aid. Please note that a hearing aid will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions. Furthermore, note that in most cases, infrequent use of a hearing aid does not permit a user to attain full benefit from it.

Consult your hearing care professional if you experience unexpected operations or events with your hearing aid.

Consult your hearing care professional if you experience unexpected operations or serious incidents with your hearing aid during use or because of its use. Your hearing care professional will support you with issue handling and, if relevant, reporting to the manufacturer and/or the national authorities.

Please note that a hearing aid will not

restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions.

Hearing aid is only a part of hearing habilitation and may need to be supplemented by auditory training and instruction in lipreading.

Furthermore, note that in most cases, infrequent use of a hearing aid does not permit a user to attain full benefit from it.

Usage of hearing aids

Hearing aids should be used only as directed and adjusted by your hearing care professional. Misuse can result in sudden and permanent hearing loss. Never allow others to wear your hearing aid, as incorrect usage could cause permanent damage to their hearing.

Choking hazards & risk of swallowing batteries and other small parts

Hearing aids, their parts, and batteries should be kept out of reach of children and

General warnings

anyone who might swallow these items or otherwise cause injury to themselves.

Batteries have occasionally been mistaken for pills. Therefore, check your medicine carefully before swallowing any pills.

If a battery or hearing aid is swallowed, see a doctor immediately and contact the National Poison Center at 1-800-222-1222 or National Battery Ingestion Hotline at 1-800-498-8666.

Battery use

Always use batteries recommended by your hearing care professional. Batteries of low quality may leak and cause bodily harm.

Never attempt to recharge your batteries, and never dispose of batteries by burning them. There is a risk that the batteries will explode.

Dysfunction

Be aware of the possibility that your hearing aid may stop working without notice. Keep this in mind when you depend on warning sounds (e.g. when you are in traffic). The hearing aids may stop functioning, for instance if the batteries have expired or if the tubing is blocked by moisture or earwax.

Active implants

The hearing aid has been thoroughly tested and characterized for human health according to international standards for human exposure (Specific Absorption Ratio - SAR), induced electromagnetic power and voltages into the human body.

The exposure values are well below international accepted safety limits for SAR, induced electromagnetic power and voltages into the human body defined in the standards for human health and coexistence with active medical

implants, such as pacemakers and heart defibrillators.

The Autophone magnet or MultiTool (which has a built-in magnet) should be kept more than 30 cm away from the implant, e.g. do not carry it in your breast pocket.

In general, please follow the guidelines recommended by the manufacturers of implantable defibrillators and pacemakers on use with magnets.

X-ray/CT/MR/PET scanning, electrotherapy and surgery

Remove your hearing aid before X-ray, CT/ MR/PET scanning, electrotherapy, surgery, etc. as your hearing aid may be damaged when exposed to strong fields.

Heat and chemicals

The hearing aid must never be exposed to extreme heat, e.g. left inside a parked car in the sun.

The hearing aid must not be dried in microwave ovens or other ovens.

The chemicals in cosmetics, hairspray, perfume, aftershave lotion, suntan lotion, and insect repellent can damage the hearing aid. Always remove your hearing aid before applying such products and allow time to dry before use.

Power instrument

Special care should be exercised in selecting, fitting and using a hearing aid whose maximum sound pressure capability exceeds 132 dB SPL (IEC 711), as there may be risk of impairing the remaining hearing of the hearing aid user.

For information on whether your hearing aid is a power instrument, see the back of this leaflet.

General warnings

Possible side effects

Hearing aids and earpieces may cause an accelerated accumulation of earwax.

The otherwise non-allergenic materials used in hearing aids may in rare cases lead to a skin irritation or other side effects.

Please seek consultation with a physician if these conditions occur.

Use on aircraft

If your hearing aid contains Bluetooth you need to activate flight mode or turn off your hearing aid on board on an aircraft, unless Bluetooth is specifically permitted by the flight personnel.

The hearing aid needs to have a pushbutton for you to be able to activate flight mode, otherwise you need to turn off your hearing aid.

Please see the model overview if your hearing aid has Bluetooth.

Connection to external equipment

The safety of the use of hearing aids when connected to external equipment (with an auxiliary input cable and/or with USB cable and/or directly), is determined by the external signal source. When connected to external equipment plugged into a wall outlet, this equipment must comply with IEC-62368-1 or equivalent safety standards.

Use of third party accessories

Only use accessories, transducers or cables supplied by the manufacturer. Incompatible accessories may result in reduced electromagnetic compatibility (EMC) of your hearing aid.

Modification of hearing aids is not allowed

Changes or modifications not expressly approved by the manufacturer will void the warranty of the equipment.

Interference

The hearing aid has been thoroughly tested for interference according to the most stringent international standards. Electromagnetic interference may occur in the vicinity of equipment with this symbol. Portable and mobile RF (radio frequency) communications equipment can affect the performance of the hearing aid. If your hearing aid is affected by electromagnetic interference, move away from the source to reduce the interference.

Troubleshooting guide

Symptom	Possible causes
No sound	Worn-out battery
	Clogged sound outlet
	Clogged microphone inlet
	Hearing aid microphone muted
Intermittent or reduced sound	Clogged sound outlet
	Moisture
	Worn-out battery
Squealing noise	Hearing aid not inserted properly
	Ear wax accumulated in ear canal
Pairing issue with Apple device	Bluetooth connection failed
	Only one hearing aid paired

Solutions

Replace the battery

Clean sound outlet* or replace the wax filter

Clean microphone inlet* or replace filter (T-Cap or O-Cap)

Unmute the hearing aid microphone

Clean sound outlet or replace the wax filter

Wipe battery and hearing aid with a dry cloth

Replace the battery

Re-insert the hearing aid

Have ear canal examined by your doctor

1) Unpair your hearing aids **2)** Turn Bluetooth off and on again. **3)** Open and close battery drawer on hearing aid. **4)** Re-pair hearing aids (visit: www.oticon.com/support).

If none of the above solutions work, consult your hearing care professional for assistance.

* According to guideline in this booklet

Cell phone

Some hearing aid users have reported a buzzing sound in their hearing aid when using cell phones, indicating that the cell phone and hearing aid may not be compatible.

The ANSI C63.19 standard determines the prediction of compatibility between a specific hearing aid and a cell phone, thus hearing aid compliance is tested according to this standard. However, demonstrating compliance according to this standard cannot guarantee that all users will be satisfied.

Whereas all hearing aids have acoustic coupling, only the larger hearing aids have the physical space for telecoil (inductive) coupling.

The hearing aid is compliant with ANSI C63.19 in both microphone and telecoil mode.

IMPORTANT NOTICE

The performance of individual hearing aids may vary with individual cell phones. Therefore, ensure you try this hearing aid with your cell phone or, if you are purchasing a new phone, be sure to try it with your hearing aid prior to purchase. For additional guidance, please ask your cell phone provider for the booklet entitled "Hearing Aid Compatibility with Digital Wireless Cell Phones."

Water & dust resistant (IP68)

Your hearing aid is dust-tight and protected against ingress of water which means it is designed to be worn in all daily life situations. Therefore you do not have to worry about sweat or getting wet in the rain. Should your hearing aid come in contact with water and stop working, please follow these guidelines:

1. Gently wipe off any water.
2. Open the battery drawer and remove the battery and gently wipe off any water in the battery drawer.
3. Let the hearing aid dry with the battery drawer left open for approximately 30 minutes.
4. Insert a new battery.

IMPORTANT NOTICE

Do not wear your hearing aid while showering or participating in water activities. Do not immerse your hearing aid in water or other liquids.

Conditions of use

Operating conditions	Temperature: +1°C to +40°C (34°F to 104°F) Relative humidity: 5% to 93%, non-condensing
Storage and transportation conditions	Temperature and humidity should not exceed the following limits for extended periods during transportation and storage: Temperature: -25°C to +60°C (-13°F to 140°F) Relative humidity: 5% to 93%, non-condensing

Warranty

Certificate

Name of owner: _____

Hearing aid professional: _____

Hearing aid professional's address: _____

Hearing aid professional's phone: _____

Purchase date: _____

Warranty period: _____ Month: _____

Model left: _____ Serial no.: _____

Model right: _____ Serial no.: _____

Battery size: _____

International warranty

Your hearing aid is covered by an international limited warranty issued by the manufacturer for a period of 12 months from the date of delivery to you. This limited warranty covers manufacturing and material defects in the hearing aid itself, but not in accessories such as batteries, tubing, speakers, earpieces and filters, etc. Problems from improper or incorrect handling or care, excessive use, accidents, repairs made by an unauthorized party, exposure to corrosive conditions, physical changes in your ear, damage due to foreign objects entering the device, or incorrect adjustments are NOT covered by the limited warranty and may void it. Under this warranty, the manufacturer will choose, in its sole discretion, whether to repair the hearing aid, or replace it with an equivalent model. The above warranty

does not affect any legal rights that you might have under applicable national legislation governing sale of consumer goods. Your hearing care professional may have issued a warranty that goes beyond the clauses of this limited warranty. Please consult him/her for further information.

If you need service

Take your hearing aid to your hearing care professional, who may be able to sort out minor problems and adjustments immediately. Your hearing care professional can assist you to obtain warranty service from the manufacturer. Your hearing care professional may charge a fee for their services.

Technical information

The hearing aid contains two radio technologies which are described below:

The hearing aid contains a radio transceiver using short range magnetic induction technology working at 3.84 MHz. The magnetic field strength of the transmitter is very weak and always below 15 nW (typically below $-40 \text{ dB}\mu\text{A}/\text{m}$ ($-12.20 \text{ dB}\mu\text{A}/\text{ft}$) at 10 meter (33 feet) distance).

The hearing aid also contains a radio transceiver using Bluetooth Low Energy (BLE) and a proprietary short range radio technology both working at ISM band 2.4 GHz. The radio transmitter is weak and always below 3 mW equal to 4.8 dBm in total radiated power.

Only use your hearing aid in areas where wireless transmission is permitted.

The hearing aid complies with international standards concerning radio transmitters, electromagnetic compatibility and human exposure. Due to the limited space available on the hearing aid relevant approval markings can be found in this booklet.

Additional information can be found in the “Technical Data sheets” on www.oticon.com

CIC only

This device contains a radio module with the following certification ID number:

FCC ID: U28-AUCIC

ITC, ITE HS & ITE FS only

This device contains a radio module with the following certification ID number:

Battery size 13:

FCC ID: U28-AUITE13

Battery size 312:

FCC ID: U28-AUITE312

Radiofrequency radiation exposure information

This device complies with FCC RF exposure limits set forth for an uncontrolled environment and has been tested for portable use.

The device must not be co-located or used in conjunction with any other antenna or transmitter.

Use of other accessories not verified by the manufacturer may not ensure compliance with FCC RF exposure guidelines.

Note: This device has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the manufacturer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The manufacturer declares that this hearing aid is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

Declaration of Conformity is available from the manufacturer.



Oticon A/S
Kongebakken 9
DK-2765 Smørum
Denmark
www.oticon.global

CE 0543



Waste from electronic equipment must be handled according to local regulations.

 **Bluetooth®**



IP68

Description of symbols used in this booklet



Warnings

Text marked with a warning symbol must be read before using the device.



Manufacturer

The device is produced by the manufacturer whose name and address are stated next to the symbol. Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.



CE mark

The device complies with Medical Device Directive 93/42/EEC. The four digits number indicates the identification of the notified body.



Electronic waste (WEEE)

Recycle hearing aids, accessories or batteries according to local regulations. Hearing aid users can also return electronic waste to their hearing care professional for disposal. Electronic equipment covered by Directive 2012/19/EU on waste and electrical equipment (WEEE).



Regulatory Compliance Mark (RCM)

The device complies with electrical safety, EMC and radio spectrum requirements for devices supplied to the Australian or New Zealand market.



IP code

Indicates the class of protections against harmful ingress of water and particulate matter according to EN 60529:1991/A1:2002. IP6X indicates total dust protection. IPX8 indicates the protection against the effects of continuous immersion in water.

**Bluetooth logo**

Registered trademark of Bluetooth SIG, Inc. where any use of such requires a license.

**Made for Apple badges**

Indicates that the device is compatible with iPhone, iPad and iPod touch.

**Radio Frequency (RF) transmitter**

Your hearing aid contains an RF transmitter.

Description of symbols used on the regulatory packaging label

**Keep dry**

Indicates a medical device that needs to be protected from moisture.

**Catalog number**

Indicates the manufacturer's catalog number so that the medical device can be identified.

**Serial number**

Indicates the manufacturer's serial number so that a specific medical device can be identified.

**Unique device identifier**

Indicates a carrier that contains unique device identifier information.

Your individual hearing aid settings

To be filled out by your hearing care professional.

Tinnitus SoundSupport: Limitation on use			
<input type="checkbox"/>	No limitation on use		
	Program	Start-up volume (Tinnitus)	Max volume (Tinnitus)
<input type="checkbox"/>	1	Max _____ hours per day	Max _____ hours per day
<input type="checkbox"/>	2	Max _____ hours per day	Max _____ hours per day
<input type="checkbox"/>	3	Max _____ hours per day	Max _____ hours per day
<input type="checkbox"/>	4	Max _____ hours per day	Max _____ hours per day

Power instrument Yes No

Settings overview for your hearing aid

Left

Right

Yes No

Volume control

Yes No

Yes No

Program shift

Yes No

Yes No

Mute

Yes No

Yes No

Tinnitus SoundSupport

Yes No

Volume control indicators

On Off

Beeps at min /max volume

On Off

On Off

Clicks when changing volume

On Off

On Off

Beeps at preferred volume

On Off

Battery indicators

On Off

Low battery warning

On Off

Summary of relevant studies

Clinical evaluations conducted by or for the manufacturer provide evidence to support the intended use and clinical benefits outlined in the IFU and demonstrate regulatory conformity. Clinical data is collected, assessed, and analyzed to support the performance of the hearing aids by validating that they provide sufficient audibility and hearing loss compensation based on best-practice prescriptive fitting rationales. The clinical data also demonstrate improved speech understanding and success with hearing aids using validated questionnaires and surveys.

Non-clinical data supporting the overall performance of the hearing aids includes software verification, electroacoustic verification, electrical and mechanical safety evaluation, electromagnetic compatibility (EMC) evaluation, and documentation of radio properties and performance. Additional information can be found in section Technical Information.

Technical Data

Measured according to American National Standard ANSI S3.22-2014
and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	108 dB SPL	108 dB SPL
	HF Average	102 dB SPL	102 dB SPL
Full-on Gain	Peak	41 dB	41 dB
	HF Average	38 dB	38 dB
Reference Test Gain		26 dB	26 dB
Frequency Range		100-9200 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	2 %	2 %
	800 Hz	2 %	2 %
	1600 Hz	2 %	2 %
Equivalent Input Noise Level	(omni)	18 dB SPL	18 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<30 ms	<30 ms

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.1 mA	1.1 mA
	Quiescent	1 mA	1 mA
Expected battery life (bat. size 10)	Hours	70-80	70-80
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	-	-
	Measured output at 10 mA/m	-	-
	Measured output at 31.6 mA/m	-	-

Technical Data

Measured according to American National Standard ANSI S3.22-2014
and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	116 dB SPL	116 dB SPL
	HF Average	113 dB SPL	113 dB SPL
Full-on Gain	Peak	47 dB	47 dB
	HF Average	46 dB	46 dB
Reference Test Gain		37 dB	37 dB
Frequency Range		100-9200 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	<2 %	<2 %
	800 Hz	2 %	2 %
	1600 Hz	<2 %	<2 %
Equivalent Input Noise Level	(omni)	18 dB SPL	18 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<45 ms	<45 ms

0 dB SPL ref. 20 μPa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.4 mA	1.4 mA
	Quiescent	1 mA	1 mA
Expected battery life (bat. size 10)	Hours	60-70	60-70
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	-	-
	Measured output at 10 mA/m	-	-
	Measured output at 31.6 mA/m	-	-

Technical Data

Measured according to American National Standard ANSI S3.22-2014
and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	109 dB SPL	109 dB SPL
	HF Average	104 dB SPL	104 dB SPL
Full-on Gain	Peak	47 dB	47 dB
	HF Average	42 dB	42 dB
Reference Test Gain		27 dB	27 dB
Frequency Range		100-7000 Hz	100-7000 Hz
Total Harmonic Distortion	500 Hz	<2 %	<2 %
	800 Hz	<2 %	<2 %
	1600 Hz	2 %	2 %
Equivalent Input Noise Level	(omni)	19 dB SPL	19 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<25 ms	<25 ms

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1 mA	1 mA
	Quiescent	1 mA	1 mA
Expected battery life (bat. size 10)	Hours	70-80	70-80
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	-	-
	Measured output at 10 mA/m	-	-
	Measured output at 31.6 mA/m	-	-

Technical Data

Measured according to American National Standard ANSI S3.22-2014
and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	118 dB SPL	118 dB SPL
	HF Average	115 dB SPL	115 dB SPL
Full-on Gain	Peak	52 dB	52 dB
	HF Average	49 dB	49 dB
Reference Test Gain		38 dB	38 dB
Frequency Range		100-9000 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	< 2 %	< 2 %
	800 Hz	< 2 %	< 2 %
	1600 Hz	2 %	2 %
Equivalent Input Noise Level	(omni)	17dB SPL	17dB SPL
Attack Time		< 5 ms	< 5 ms
Release Time		< 35 ms	< 35 ms

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.3 mA	1.3 mA
	Quiescent	1 mA	1 mA
Expected battery life (bat. size 10)	Hours	60-70	60-70
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	-	-
	Measured output at 10 mA/m	-	-
	Measured output at 31.6 mA/m	-	-

Technical Data

ITC/ITE HS/ITE FS

75

Measured according to American National Standard ANSI S3.22-2014

and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	108 dB SPL	108 dB SPL
	HF Average	103 dB SPL	103 dB SPL
Full-on Gain	Peak	45 dB	45 dB
	HF Average	41 dB	41 dB
Reference Test Gain		27 dB	27 dB
Frequency Range		100-7500 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	<2 %	<2 %
	800 Hz	<2 %	<2 %
	1600 Hz	<2 %	<2 %
Equivalent Input Noise Level	(omni/dir)	16/27 dB SPL	16/27 dB SPL
HF Average SPLITS	(left/right ear)	83/83 dB SPL	83/83 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<25 ms	<25 ms

0 dB SPL ref. 20 μPa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.8 mA	1.8 mA
	Quiescent	1.7 mA	1.7 mA
Expected battery life (bat. size 312/13)	Hours	55-60/100-115	55-60/100-115
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	69 dB SPL	69 dB SPL
	Measured output at 10 mA/m	89 dB SPL	89 dB SPL
	Measured output at 31.6 mA/m	100 dB SPL	100 dB SPL

Technical Data

ITC/ITE HS/ITE FS

85

Measured according to American National Standard ANSI S3.22-2014

and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	116 dB SPL	116 dB SPL
	HF Average	112 dB SPL	112 dB SPL
Full-on Gain	Peak	50 dB	50 dB
	HF Average	47 dB	47 dB
Reference Test Gain		35 dB	35 dB
Frequency Range		100-8800 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	< 2 %	< 2 %
	800 Hz	< 2 %	< 2 %
	1600 Hz	< 2 %	< 2 %
Equivalent Input Noise Level	(omni/dir)	15/27 dB SPL	15/27 dB SPL
HF Average SPLITS	(left/right ear)	92/92 dB SPL	92/92 dB SPL
Attack Time		< 5 ms	< 5 ms
Release Time		< 35 ms	< 35 ms

0 dB SPL ref. 20 μPa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.9 mA	1.9 mA
	Quiescent	1.7 mA	1.7 mA
Expected battery life (bat. size 312/13)	Hours	50-60/95-115	50-60/95-115
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	75 dB SPL	75 dB SPL
	Measured output at 10 mA/m	95 dB SPL	95 dB SPL
	Measured output at 31.6 mA/m	106 dB SPL	106 dB SPL

Technical Data

ITC/ITE HS/ITE FS

90

Measured according to American National Standard ANSI S3.22-2014

and ANSI S3.55-2014/Part 5

Supply voltage: Zinc Air 1.4 Volt

0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	120 dB SPL	120 dB SPL
	HF Average	116 dB SPL	116 dB SPL
Full-on Gain	Peak	55 dB	55 dB
	HF Average	50 dB	50 dB
Reference Test Gain		39 dB	39 dB
Frequency Range		100-7900 Hz	100-7500 Hz
Total Harmonic Distortion	500 Hz	<2 %	<2 %
	800 Hz	<2 %	<2 %
	1600 Hz	<2 %	<2 %
Equivalent Input Noise Level	(omni/dir)	15/27 dB SPL	15/27 dB SPL
HF Average SPLITS	(left/right ear)	96/96 dB SPL	96/96 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<25 ms	<25 ms

0 dB SPL ref. 20 μPa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.8 mA	1.8 mA
	Quiescent	1.7 mA	1.7 mA
Expected battery life (bat. size 312/13)	Hours	55-60/105-115	55-60/105-115
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	79 dB SPL	79 dB SPL
	Measured output at 10 mA/m	99 dB SPL	99 dB SPL
	Measured output at 31.6 mA/m	109 dB SPL	109 dB SPL

Technical Data

ITC/ITE HS/ITE FS **100**

Measured according to American National Standard ANSI S3.22-2014
and ANSI S3.55-2014/Part 5
Supply voltage: Zinc Air 1.4 Volt



0 dB SPL ref. 20 μ Pa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
OSPL90	Peak	125 dB SPL	125 dB SPL
	HF Average	122 dB SPL	122 dB SPL
Full-on Gain	Peak	63 dB	63 dB
	HF Average	58 dB	58 dB
Reference Test Gain		45 dB	45 dB
Frequency Range		100-7100 Hz	100-7100 Hz
Total Harmonic Distortion	500 Hz	<2 %	<2 %
	800 Hz	<2 %	<2 %
	1600 Hz	<2 %	<2 %
Equivalent Input Noise Level	(omni/dir)	15/28 dB SPL	15/28 dB SPL
HF Average SPLITS	(left/right ear)	103/103 dB SPL	103/103 dB SPL
Attack Time		<5 ms	<5 ms
Release Time		<15 ms	<15 ms

0 dB SPL ref. 20 μPa		Oticon Opn 1	Oticon Opn 2 & 3 Oticon Siya 1 & 2
Battery Consumption	Typical	1.8 mA	1.8 mA
	Quiescent	1.7 mA	1.7 mA
Expected battery life (bat. size 312/13)	Hours	50-60/90-115	50-60/90-115
Latency		8.2 ms	8.2 ms
Maximum Induction Coil Sensitivity	Measured output at 1 mA/m	89 dB SPL	89 dB SPL
	Measured output at 10 mA/m	109 dB SPL	109 dB SPL
	Measured output at 31.6 mA/m	119 dB SPL	119 dB SPL



Oticon Opn 1

Oticon Opn 2 & 3 Oticon Siya 1 & 2



OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0



Frequency response
Input: 60 dB SPL
Technical setting: NO





Input-Output
Technical Setting: NO





Oticon Opn 1

Oticon Opn 2 & 3 Oticon Siya 1 & 2



OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0



Frequency response
Input: 60 dB SPL
Technical setting: N0





Input-Output
Technical Setting: N0





Oticon Opn 1

Oticon Opn 2 & 3 Oticon Siya 1 & 2



OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0



Frequency response
Input: 60 dB SPL
Technical setting: NO





Input-Output
Technical Setting: NO





Oticon Opn 1

Oticon Opn 2 & 3 Oticon Siya 1 & 2



OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0



Frequency response
Input: 60 dB SPL
Technical setting: N0





Input-Output
Technical Setting: N0



Oticon Opn 1



Oticon Opn 2 & 3 Oticon Siya 1 & 2

OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0





Frequency response
Input: 60 dB SPL
Technical setting: NO

— Acoustic input
- - - Magnetic input





Input-Output
Technical Setting: NO



Oticon Opn 1



Oticon Opn 2 & 3 Oticon Siya 1 & 2

OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0





Frequency response
Input: 60 dB SPL
Technical setting: N0

— Acoustic input
- - - Magnetic input





Input-Output
Technical Setting: N0



Oticon Opn 1



Oticon Opn 2 & 3
Oticon Siya 1 & 2

OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0





Frequency response
Input: 60 dB SPL
Technical setting: NO

— Acoustic input
- - - Magnetic input



Input-Output
Technical Setting: NO





266682US/2023.02.07 / v1

Oticon Opn 1



Oticon Opn 2 & 3 Oticon Siya 1 & 2

OSPL90 - Output Sound Pressure Level
Input: 90 dB SPL
Technical setting: A0





Frequency response
Input: 60 dB SPL
Technical setting: N0

— Acoustic input
- - - Magnetic input



Input-Output
Technical Setting: N0





oticon
life-changing **technology**