


MedRx®

www.medrx-diagnostics.com



8 user-friendly modules. Clean, modern design. Simple to use and easy to customize for your clinic.

The New MedRx Studio Software

One Great Software for all the MedRx Products

MedRx®

888-392-1234

To Schedule An Online Demo




The New MedRx Studio Software

8 user-friendly modules


Incorporate the New
MedRx Studio Software
into your daily practice to
deliver quality professional
hearing healthcare.

This clean, modern
design is easy to use and
functions as your complete
Audiological suite.




DUAL CHANNEL AUDIOMETRY

Perform air conduction, bone conduction, masking and speech testing using the mouse and/or keyboard with ease.



REAL EAR MEASUREMENTS


IEC/ANSI standard REM functions are provided for: REUR, REUG, REIG, REAG, REAR, REOR, REOG and RECD.



HEARING INSTRUMENT TESTING

The IEC/ANSI standard tests are available in an automatic sequence or for select measurements.


LSM



LIVE SPEECH MAPPING

Live Speech Mapping (LSM) a quick and successful approach to hearing instrument fitting.


IMP



TYMPANOMETRY

This function provides clear graphic and numeric displays of tympanogram and acoustic reflex data.

HLS
MHA



HLS/MHA

Easily simulate hearing loss and hearing instruments for patients.


VO



VIDEO OTOSCOPE

Capture, view and store digital images. Retrieve results easily for follow-up visits.

VRA



ivRA

A new form of visual reinforcement audiometry. A2D+ Audiometer, ARC and AWRC compatible.



STUDIO

Software

The MedRx Studio Software is easy to use and equips you with excellent counseling tools.



MedRx[®]

Good Things Come In Small Packages

MedRx, Inc

1200 Starkey Rd. #105

Largo, FL 33771 U.S.A.

888-392-1234 • 727-584-9600

Fax: 727-584-9602

Email: sales@medrx-diagnostics.com

www.medrx-diagnostics.com

MedRx Minimum Computer Specs:

Windows[®] PC compatible computer, Intel[™] i5, 2.0 GHz or better. 4 GB RAM. 20 GB free hard drive space.

Available 2.0 USB Port. Windows 10 or 11 Professional (32 Or 64-Bit).

MX-STUDIO-DS-1 Rev. 4

