

OTICON | CROS

The Oticon CROS family

Help your patients stay in the conversation

life-changing
technology

Oticon CROS and CROS PX transmitters are the perfect transmitters to give your patients with single-sided deafness confidence and peace of mind throughout their day. With advanced TwinLink™ dual-streaming technology, the Oticon CROS family offers patients access to all relevant sounds in challenging listening environments while also providing the ability to stream sound simultaneously.




Oticon CROS

Oticon CROS PX

The Oticon CROS Family

- Offers patients the choice of a rechargeable (Oticon CROS PX) or non-rechargeable (Oticon CROS) transmitter to match their personal preference and hearing aid style.
- Both Oticon CROS and Oticon CROS PX offer forwards and backwards compatibility with Oticon More™, Oticon Zircon, Oticon Play PX, Opn S™*, Ruby, Opn Play™*, Xceed and Xceed Play.
- Includes TwinLink, an innovative technology featuring simultaneous Near-Field Magnetic Induction (NFMI) and 2.4 GHz Bluetooth® Low Energy streaming, making it possible to carry on a conversation while streaming sound.
- Can be used as either a CROS or BICROS fitting.
- Compatible with miniRITE, miniRITE T, miniRITE R, miniBTE T, miniBTE R, BTE, BTE PP, BTE SP and BTE UP.
- Both styles available in the 7 standard colors with Hear Pink also available for Oticon CROS PX.

*Except Opn S 3 and Opn Play 2





Oticon CROS family compatibility



Oticon CROS
miniRITE T



Oticon CROS PX
miniRITE R

Oticon More 1*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon More 2*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon More 3*: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon Opn S 1**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn S 2**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn S 3**: miniRITE, miniRITE T, miniRITE R, BTE PP

Oticon Xceed 1**: BTE SP, BTE UP

•

Oticon Xceed 2**: BTE SP, BTE UP

•

Oticon Xceed 3**: BTE SP, BTE UP

•

NEW **Oticon Zircon 1***: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

NEW **Oticon Zircon 2***: miniRITE R, miniRITE T, miniBTE R, miniBTE T

•

Oticon Ruby 1***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

•

Oticon Ruby 2***: miniRITE, miniRITE T, miniRITE R, BTE, BTE PP

•

NEW **Oticon Play PX 1***: miniRITE T, miniRITE R, miniBTE R, miniBTE T

•

NEW **Oticon Play PX 2***: miniRITE T, miniRITE R, miniBTE R, miniBTE T

•

Oticon Opn Play 1**: miniRITE, miniRITE T, miniRITE R, BTE PP

•

Oticon Opn Play 2**: miniRITE, miniRITE T, miniRITE R, BTE PP

Oticon Xceed Play 1**: BTE SP, BTE UP

•

Oticon Xceed Play 2**: BTE SP, BTE UP

•

• Indicates feature included

* Firmware 1.1.0 and onwards

** Firmware 8.0 and onwards

*** Firmware 1.0 and onwards

21-176650 15500-0567/12.21



oticon.com

Oticon is part of the Demant Group.

oticon
life-changing technology