

Feature overview

	Own SI 1	Own SI 2	Own SI 3	Own SI 4
Speech understanding & listening ease				
MoreSound Intelligence™ 3.0	Level 1	Level 2	Level 3	Level 4
Environment classifier	5 configurations	5 configurations	3 configurations	Not adjustable
Neural Noise Suppression Difficult/Easy	12 dB / 6 dB	10 dB / 4 dB	8 dB / 2 dB	6 dB / 0 dB
Sound Enhancer	3 configurations	2 configurations	1 configuration	1 configuration
MoreSound Amplifier™ 3.0	•	•	•	•
SuddenSound Stabilizer	6 configurations	5 configurations	4 configurations	2 configurations
MoreSound Optimizer™	•	•	•	•
Feedback shield	•	•	•	•
Spatial Sound™ ¹	o	o	o	
Soft Speech Booster	•	•	•	•
Frequency lowering, Speech Rescue™	•	•	•	•
Sound quality				
Clear Dynamics	•	•		
Better-Ear Priority ¹	o	o	o	
Fitting Bandwidth ²	10 kHz	8 kHz	8 kHz	8 kHz
Processing Channels	64	48	48	48
Personalization & optimized fitting				
Fitting Bands	24	20	18	14
Adaptation Management	•	•	•	•
Fitting formulas	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5	VAC+, NAL-NL1/ NAL-NL2, DSL v5
Audible Contrast Threshold (ACT™) prescription	•	•	•	•
Tinnitus SoundSupport™ ³	o	o	o	o

¹Requires NFMI

²Bandwidth accessible for gain adjustments during fitting


³Requires NFMI and push-button

• Default

o Optional features only available for CIC

OTICON | Own SI

Product comparison




oticon
life-changing technology

Oticon Own SI™ styles and fitting options


Oticon Own SI IIC and CIC are our smallest in-the-ear styles. The hearing aids are powered by the Sirius™ platform and include Oticon BrainHearing™ technology. IIC and CIC are discreet hearing aids with the IIC being invisible in most ears. Both styles use disposable batteries. Wireless accessories are not compatible with Own SI IIC and CIC.



Oticon Own SI 1




75




90


Oticon Own SI 2 | 3 | 4



75



90



C0001
Beige



C0002
Light Brown



C0003
Medium Brown



C0004
Dark Brown



C0005
Black





Transparent



Transparent
Blue



Transparent
Red

Style	Battery size	Fitting level	NMFI (Binaural Coordination Only)	Microphones	Push button	IP68 - Water and dust resistant (Hardware certification)
 IIC	10	75 90		1		✓
 CIC	10	75 90	o	1	o	✓

✓ Default

o Optional