

Your Information	Step 1: Account Information					Step 2: Patient Information								
	Ship to Account:					First Name:								
	Account Name:					Last Name:								
	Address:					Age:								
						Last 4 digits of Social Security #:								
	City:		State:		Zip:		Audiometric Information (required)							
	Phone #:					Hz:	250	500	1k	2k	3k	4k	6k	8k
	Contact Name:					Right:								
	Email (required):					Left:								
	Bill to Account:					Step 3: Fitter's Information								
Date:			Fitting Date:		Fitter's Name:									
Purchase Order #:					Fitter's Email:									

Existing Instruments	
Please indicate existing instrument information:	Please indicate existing instrument information:
Model: _____ Serial #: _____	Model: _____ Serial #: _____
What type of speaker are you fitting?	What type of speaker are you fitting?
<input type="checkbox"/> miniFit 60 <input type="checkbox"/> miniFit 85 <input type="checkbox"/> miniFit 100 <input type="checkbox"/> miniFit 105	<input type="checkbox"/> miniFit 60 <input type="checkbox"/> miniFit 85 <input type="checkbox"/> miniFit 100 <input type="checkbox"/> miniFit 105

BTE Mold Order													
Type	L	R	Style	L	R	Options	L	R	Slim Tube	L	R	0.9mm†	1.3mm†
BTE Mold	<input type="checkbox"/>	<input type="checkbox"/>	Canal	<input type="checkbox"/>	<input type="checkbox"/>	Name/initials on mold:	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size -1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Slim Tube Mold (Corda)	<input type="checkbox"/>	<input type="checkbox"/>	Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	_____	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size 0	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Material</b>	<b>L</b>	<b>R</b>	Half Shell	<input type="checkbox"/>	<input type="checkbox"/>	Print L&R on mold	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size 1	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Hard (Clear STD)	<input type="checkbox"/>	<input type="checkbox"/>	Full Shell	<input type="checkbox"/>	<input type="checkbox"/>	Make Canal Tip red & blue (Hard Only)	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size 2	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Soft silicone (Clear STD)	<input type="checkbox"/>	<input type="checkbox"/>	Helix Lock (n/a in silicone)	<input type="checkbox"/>	<input type="checkbox"/>	Blue/red/white dot on mold <sup>1</sup>	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size 3	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
OtoTherm™ (Clear STD)	<input type="checkbox"/>	<input type="checkbox"/>	Skeleton (n/a in silicone)	<input type="checkbox"/>	<input type="checkbox"/>	Removal String	<input type="checkbox"/>	<input type="checkbox"/>	Corda miniFit tube size 4	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
			Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	Removal String heavy	<input type="checkbox"/>	<input type="checkbox"/>	<b> tubing</b>	<b>L</b>	<b>R</b>		
			Semi Skeleton (n/a in silicone)	<input type="checkbox"/>	<input type="checkbox"/>	<b>Venting</b>	<b>L</b>	<b>R</b>	13 Medium tube	<input type="checkbox"/>	<input type="checkbox"/>		
			<b>Length</b>	<b>L</b>	<b>R</b>	Small (0.8)	<input type="checkbox"/>	<input type="checkbox"/>	13 DRY Medium tube	<input type="checkbox"/>	<input type="checkbox"/>		
			Short Canal	<input type="checkbox"/>	<input type="checkbox"/>	Medium (1.4)	<input type="checkbox"/>	<input type="checkbox"/>	13 Heavy tube	<input type="checkbox"/>	<input type="checkbox"/>		
			Medium Canal	<input type="checkbox"/>	<input type="checkbox"/>	Large (2.4)	<input type="checkbox"/>	<input type="checkbox"/>	13 DRY Heavy tube	<input type="checkbox"/>	<input type="checkbox"/>		
			Long Canal	<input type="checkbox"/>	<input type="checkbox"/>	Max	<input type="checkbox"/>	<input type="checkbox"/>	13 Double-wall tube	<input type="checkbox"/>	<input type="checkbox"/>		
			As Marked	<input type="checkbox"/>	<input type="checkbox"/>	SAV	<input type="checkbox"/>	<input type="checkbox"/>					
						Pressure	<input type="checkbox"/>	<input type="checkbox"/>					
						No Vent	<input type="checkbox"/>	<input type="checkbox"/>					

<sup>1</sup>N/A Soft  
<sup>†</sup>Corda thin tubes 0.9mm and 1.3mm are only compatible with Oticon Real miniBTEs, Oticon More miniBTEs and Opn S Power Plus BTEs.

Domes for Corda Thin Tubes										Notes	
	Open*		OpenBass**		Double Vent Bass (1.4)		Single Vent Bass (0.7)*		Power		<input type="checkbox"/> Order contains <b>contaminated impression</b> . Handlers use gloves.
	L	R	L	R	L	R	L	R	L	R	
5mm	N/A	N/A	<input type="checkbox"/>	<input type="checkbox"/>	N/A	N/A	N/A	N/A	N/A	N/A	Power Molds default to 1.4mm vents unless otherwise requested. Impression will be destroyed unless otherwise requested. <b>Special Instructions:</b> _____ _____ _____
6mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
8mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
10mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
12mm	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

\* Compatible only with Opn S 1 Plus Power  
 \*\* Compatible only for use with Oticon Real™ and Oticon More™ with the exception of 5mm OpenBass which is equivalent to the 5mm Open dome.

Standard (Speakers 60, 85)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™	Venting	Removal String		
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Please choose Venting selection options below.	STD	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>			<input type="checkbox"/>	
	Hollow Mold (LiteTip style)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)		OtoTherm™	Venting	Removal String
Hollow	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>	STD			
Hollow with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>	<input type="checkbox"/>			
Hollow with Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	<input type="checkbox"/>	<input type="checkbox"/>			
Embedded									
MicroShell for 60 & 85 Speakers (embedded speaker)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™	Venting	Options		
Canal	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-		MicroShell (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 60 <input type="checkbox"/> 85	<input type="checkbox"/> R <input type="checkbox"/> 60 <input type="checkbox"/> 85
Canal with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-			<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5	<input type="checkbox"/> R <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
Canal with Helix Lock	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-		Receiver Length	<input type="checkbox"/> Clear <input type="checkbox"/> Beige	<input type="checkbox"/> Black <input type="checkbox"/> Lt. Brown
Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-			<input type="checkbox"/> Med. Brown <input type="checkbox"/> Dark Brown	<input type="checkbox"/> Blue <input type="checkbox"/> Red
Semi Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-		Shell and Faceplate color		
Half Skeleton	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	-Not Available-	-Not Available-				
Power Receiver Mold (100 and 105 speakers)	L	R	Hard material acrylic (clear only)	Soft material silicone (no tint)	OtoTherm™	Venting	Options		
Power Receiver Mold	<input type="checkbox"/>	<input type="checkbox"/>	STD	-Not Available-	-Not Available-		Power Receiver Mold (acrylic only)	<input type="checkbox"/> L <input type="checkbox"/> 100 <input type="checkbox"/> 105	<input type="checkbox"/> R <input type="checkbox"/> 100 <input type="checkbox"/> 105
Power Receiver Mold with Canal Lock	<input type="checkbox"/>	<input type="checkbox"/>	STD	-Not Available-	-Not Available-			Power Receiver Wire Length	<input type="checkbox"/> L <input type="checkbox"/> 1 <input type="checkbox"/> 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/> 5
							Shell and Cover Plate Color		<input type="checkbox"/> Clear <input type="checkbox"/> Beige <input type="checkbox"/> Lt. Brown <input type="checkbox"/> Med. Brown <input type="checkbox"/> Dark Brown <input type="checkbox"/> Black
Removal String		<input type="checkbox"/> Yes <input type="checkbox"/> No							

Venting Selection															
All RITE molds except Hollow (LiteTip) Mold								Hollow (LiteTip) Mold RITE							
Venting based on audio-canal size <input type="checkbox"/>								Venting based on audio-canal size <input type="checkbox"/>							
Vent Size	Small (0.8mm)	Medium (1.4mm)	Large (2.4mm)	Max**	Pressure	SAV	No vent	Vent Size	Small (0.7mm long)	Medium (0.7mm short)	Large (1.2mm)	Max**	Pressure	SAV	No vent
L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	L	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

\*\*As large as possible.

**Retention Locks**

Available options vary by style and material selected.

Canal Lock	Helix Lock	Skeleton Mold	Half Skeleton	Semi Skeleton
	(Not available in soft material) 			(Not available in soft material) 

**Wax Protection System**

Custom molds (except OtoTherm™, MicroShell, and Soft Silicone MicroMold) use the ProWax system. All miniFit receivers use the ProWax miniFit wax system.



\*OtoTherm™ is flexible material that is rigid (Shore A 90) at room temperature and flexible (Shore A 70) at body temperature. These characteristics allow for ease of insertion, increased comfort, and optimal sealing.

Please note: This material requires gentle warming/softening of the mold with a hair dryer prior to insertion or removal of the receiver.