

Oticon Government Services

life-changing  
technology

# Selection Guide Book

Check this out!




**NEW**

Scan the QR code above to add the Oticon Patient Support Hotline to your contacts.



Phone compatibility check

May 2025



**oticon**  
life-changing technology

## RITE Styles

### Oticon Intent 1 miniRITE R (lithium-ion)

- Always-on DNN 2.0
- 4D Sensor technology
- Bluetooth® LE Audio
- SuddenSound Stabilizer

### Oticon Real 1 miniRITE R (lithium-ion)

- Always-on DNN
- SuddenSound Stabilizer
- Toggle Button

### Oticon Play PX 1 miniRITE R (lithium-ion)

#### Oticon Opn S 1 miniRITE R (lithium-ion)

Meets SCIF requirements

### Oticon Real 1 miniRITE T (312 zinc-air)

- Always-on DNN
- SuddenSound Stabilizer
- Toggle Button

60 85 100 105

60 85 100 105

60 85 100 105

60 85 100 105

### Oticon Real 1 miniBTE R (lithium-ion)

- Always-on DNN
- SuddenSound Stabilizer
- BTE or Corda tubing

### Oticon Real 1 miniBTE T (312 zinc-air)

- Always-on DNN
- SuddenSound Stabilizer
- BTE or Corda tubing

### Opn S 1 BTE Plus Power 13

Meets SCIF requirements

### Xceed 1 BTE Super Power 13

Meets SCIF requirements

### Xceed 1 BTE Ultra Power 675

Meets SCIF requirements

85

85


105

110

120




## Intent Fitting Ranges




60

- OpenBass dome
- Mold, Bass & Power dome




85



100


- MicroShell Detect 100, Bass & Power dome



105


- MicroShell Detect 105

## RITE Fitting Ranges




60

- Open Bass/Open dome
- Mold, Bass & Power dome




85



100


- Bass, Power dome & Power Receiver Mold



105


- Power Receiver Mold

## BTE Fitting Ranges




85

- Corda miniFit standard thin tube




105

- Corda miniFit Power thin tube



110

- Hook



120

## Oticon RITE and BTE Colors

C044

Silver

C092

Steel Grey

C091

Silver Grey

C090

Chroma Beige

C116\*\*

Honey Beige

C094

Terracotta

C093

Chestnut Brown

C063

Diamond Black

C079

Hear Pink\*

C115\*\*

Sky Blue

C114†

Olive Green

## Oticon Play PX Colors Only

C048

Emerald Green

C058

Aquamarine

C047

Cool Blue

C045

Purple

C046

Red

C057

Power Pink

C079

Baby Pink

\*Only available for Oticon Real and CROS PX.

\*\*Only available for Oticon Intent.

†Only available for Oticon Real miniRITE R.

## Custom Styles

<b>Oticon Own 1 IIC</b> Meets SCIF requirements	<b>Oticon Own 1 CIC</b> Meets SCIF requirements	<b>Oticon Own 1 ITC</b> Meets SCIF requirements if ordered without 2.4GHz	<b>Oticon Own 1 ITE HS</b> Meets SCIF requirements if ordered without 2.4GHz	<b>Oticon Own 1 ITE FS</b> Meets SCIF requirements if ordered without 2.4GHz
--	--	--	---	---

75 90

75 90

75 90 100

75 90 100

75 90 100



## Oticon CROS Family

**TwinLink Dual Streaming**  
 Backward and forward compatibility • Two battery options

### CROS Compatibility

Oticon Intent 1	Oticon Real 1	Oticon More 1	Opn S 1
		Oticon Play PX 1	Xceed 1





Oticon CROS



Oticon CROS PX

## Custom Fitting Ranges

IIC                      CIC                      ITC                      ITE HS                      ITE FS



## Chargers

### Oticon Intent™ Chargers

- Quick charge in 15 minutes
- Advanced contact charging
- Fitting during charging available
- Power Bank<sup>1</sup>



Oticon Intent SmartCharger<sup>2</sup>  
268271



Oticon Intent Charger<sup>2</sup>  
273377

### Oticon SmartChargers for Oticon Real™, Oticon More™ and Oticon Play PX<sup>3</sup>

- Protective Lid
- Power Bank



Oticon miniRITE SmartCharger  
225678



Oticon miniBTE SmartCharger  
226282

### Oticon Desktop Chargers for Oticon Real, More, Play PX and Opn S™

- Easy insertion and removal
- Contact Free<sup>4</sup>
- Quick charge<sup>4</sup>



Oticon miniRITE Charger  
200294



Oticon miniBTE Charger<sup>5</sup>  
226352

## Oticon Custom Colors

**C001**  
Beige

**C002**  
Light Brown

**C003**  
Medium Brown

**C004**  
Dark Brown

**C005**  
Black

Transparent\*

Red\*

Blue\*

\*Not available as a faceplate color.

<sup>1</sup>Oticon Intent SmartCharger Only.

<sup>2</sup>Oticon Intent Charger only compatible with Oticon Intent.


<sup>3</sup>Oticon miniRITE SmartCharger is compatible with Oticon More devices on Firmware 1.1.1 or newer.

Oticon miniRITE SmartCharger is not compatible with Opn S miniRITE R.

<sup>4</sup>Features for both Oticon SmartCharger and Desktop charger.

<sup>5</sup>Not compatible with Opn S.

## miniFit Earpiece Compatibility Overview




\*Opn S only.

\*\*Not compatible with Opn S except for 5mm OpenBass dome.

### RITE molds

### BTE Molds




Now available with  
Extended Receiver Tubes!

\*Compatible with miniFit receivers 60 and 85.

\*\*Oticon Intent only. Compatible with Oticon Intent miniFit Detect receivers 60, 85, 100 and 105.

\*\*\*Not available for Oticon Intent.

### Standard Earmold Styles\*



\*Reference pages 38-39 of pocket guide for more information on these earmold styles.

### Special Earmold Styles\*



\*Reference pages 38-39 of pocket guide for more information on these earmold styles.

### Wireless Accessories\*



#### EduMic

- Automatic pairing
- 3.5mm input jack
- Telecoil
- Easy-to-use remote microphone



#### TV Adapter 3.0

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 50 ft
- Digital/RCA connection




#### ConnectClip

- Make and receive hands-free calls to both iPhone® and Android™ devices
- Remote Microphone
- Control the volume and change programs



#### Remote Control 3.0

- Adjust volume
- Switch programs
- Mute



#### Phone Adapter 2.0

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 65 ft
- ConnectClip required



#### BTDB800a Bluetooth Dongle

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required

Apple, the Apple logo, iPhone, iPad, Vision Pro, and Mac are trademarks of Apple Inc., registered in the U.S. and other countries. Android™ is a trademark by Google LLC. The Easy LE Adapter is manufactured by Minami Acoustics Limited. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. Any use of such marks by Demant and Minami Acoustics Limited is under license. Other trademarks and trade names are those of their respective owners.

\*Only available in styles with 2.4 GHz.