

Oticon Real™ hearing aids



Experience the real sounds of life

Oticon Real with BrainHearing™ technology works more like the brain because it learned through experience. BrainHearing technology gives your brain access to virtually all real-world sounds.



Oticon Real:

- Offers technology that supports how the brain naturally processes sounds for better access to the full sound scene.
- Provides better access to speech and exceptional speech clarity, so you can follow conversations with more awareness, focus and control.
- Significantly reduces listening effort* with the Wind & Handling Stabilizer. Proven to increase speech clarity in windy situations** and minimize distracting handling noises (i.e. your hair or glasses touching the hearing aid).

Connect to your digital world

Staying connected has never been easier with Oticon Real.

- Enjoy the convenience of two-way hands-free*** communication with select iPhone® and iPad® devices and direct streaming from select Android™ devices****
- Wirelessly stream someone's voice to your hearing aids with ConnectClip and EduMic remote microphones
- Stream sound directly from your TV with the Oticon TV adapter
- Use the new Oticon Companion app to easily control and personalize your hearing aids

Real power. Real styles.

Oticon Real has a distinct, elegant design and is available in 4 styles, 9 colors† and 2 convenient charger options that make it easy to stay powered all day.††

- miniRITE R & miniBTE R styles are available for those who want the convenience of a rechargeable hearing aid
- miniRITE T & miniBTE T styles are ideal for those who prefer disposable batteries
- Introducing new olive green color, available in miniRITE R style
- Choice of a sleek desktop charger or the portable SmartCharger for on-the-go rechargeability



Apple, the Apple logo, iPhone®, iPad®, and iPod touch® are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android™, Google Play, and the Google Play logo are trademarks of Google LLC.

*Santurette et al. (2023) SuddenSound Stabilizer – Evidence and user benefits. Oticon Whitepaper

**Gade et al. (2023) Wind & Handling Stabilizer – Evidence and user benefits. Oticon Whitepaper

***Hands-free communication is available with iPhone® 11 or later running iOS 15.2 or later, and iPad® running iPadOS® 15.2 or later.

****Android™ devices need to support Audio Streaming for Hearing Aids (ASHA) to allow direct streaming to Oticon Real. Please visit oticon.com/support/compatibility for more information.

†Olive green is only available in miniRITE R style.

††Lithium-ion battery performance varies depending on hearing loss, lifestyle and streaming behavior.



oticon.com

Oticon is part of the Demant Group.

oticon
life-changing technology