

HearLink 9050

Premium digital comparison chart

When purchasing hearing aids, it is important to compare features, service, and price in your buying decision. Use this chart to help you assess Philips HearLink against other hearing aids and technology available.

HearLink 9050 features	HearLink 9050 specifications	Competitive specifications
SoundMap 3 Amplification		
Frequency bandwidth	80 Hz to 10 kHz	
Channels	64	
Extended dynamic range	●	
Low frequency enhancement	●	
Frequency Lowering	●	
SoundMap Noise Control		
Speech Clarifier	3 Options	
Comfort Control	4 Options	
SoundMap Noise Control Transition	4 Options	
Directionality		
Dynamic Directionality	3 Options	
Fixed Directionality	●	
Adaptive Directionality	●	
Pinna Mode	2 Options	
AI Noise Reduction Mode		
Noise Reduction Mode	5 Options	
Specific Noise Management		
SoundGuide motion sensor Spatial Awareness	●	
Soft Noise Management	●	
SoundProtect Transient Noise Reduction	6 Options	
SoundProtect Wind Noise Management	●	
Binaural Noise Management	●	
SoundMap Feedback Canceller		
Strength Control	●	
SoundTie 3 Connectivity: Hearing Aid*		
NFMI and 2.4 GHz Bluetooth® Low Energy	●	
LE Audio Auracast™ broadcast audio ready	●	
iOS® and Android™ direct streaming	●	
Hands-free comm. for iPhone®, iPad® and Android™	●	
Push button tap control functionality	●	
Philips wireless accessories	●	
SoundTie 3 Connectivity: App*		
HearLink 2 app: Remote Control + Remote Fitting	●	
Hearing Journal	●	
Sound Equalizer	●	
Remember volume settings	●	
Apple Watch® compatibility	●	
Binaural coordination		
Binaural Volume and Program Change	●	
Programming options		
Fitting bands	24	
Environments	13	
Manual listening programs	4	
HiFi Music program	●	
Airplane program	●	
Data Logging	●	
Connection Count	●	
Adaptation Manager	●	
Remote Fitting	●	
Notify Me	●	
Additional items		
Rechargeability	Yes – miniRITE	
Service warranty	Hearing aid: 3-year warranty; Remote Control: 2-year warranty	
Price \$		Price \$

© 2024 Demant. Philips and the Philips Shield Emblem are registered trademarks of Koninklijke Philips N.V. iPhone and iPad are trademarks of Apple Inc., registered in the U.S. and other countries. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Demant is under license. Other trademarks and trade names are those of their respective owners.